



## Anti Galanin (20-29)(Rat) Serum

Cat. No. YII-Y186-EX      Lot No. 383160701

**Description:** This antiserum was raised in a rabbit by immunization with a carrier free synthetic preprogalanin (89-124)(rat) peptide. The product vial contains 50 $\mu$ L of the titled antiserum obtained by lyophilizing its 0.001M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with reprogalanin (rat).

**Immunogen:** Synthetic galanin (20-29) (rat), carrier free

**Amino Acid Sequence of galanin (20-29)(rat)<sup>1)</sup>:**

**Product Form:** Lyophilized unpurified serum **Size:** 50  $\mu$ L

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN<sub>3</sub> 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstituted, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

**Suggested Working Dilution Range:** 1:3,000 (final dilution ~1:21,000) for radioimmunoassay; 1:500-2,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

**Specificity** (based on radioimmunoassay): Preproalanin (rat) 100%, preprogalanin (human) 80%, galanin (rat) 0%

**Positive Control** (immunohistochemistry): Rat pyroric spincter

**Species Tested:** Rat

### REFERENCES:

1) L.M. Kaplan, E.R. Spindel et al., Tissue-specific expression of the rat galanin gene. Proceedings of National Academy of Sciences, USA 85: 1065-1069, 1988

### FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

#### Distributor



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

#### Manufacturer



2480-1 Awakura, Fujinomiya-shi,  
Shizuoka 418-0011, JAPAN  
Phone: +81-544-22-2771  
Fax: +81-544-22-2770  
URL : [www.yanaihara.co.jp](http://www.yanaihara.co.jp)  
E-mail : [ask@yanaihara.co.jp](mailto:ask@yanaihara.co.jp)