



## Anti PTHrP (15-34)-NH<sub>2</sub> Serum

Cat. No. YII-Y202-EX      Lot No. 092271024

**Description:** This antiserum was raised in a rabbit by immunization with a carrier free PTHrP (15-34)-NH<sub>2</sub> peptide. The product vial contains 50 µL of the titled antiserum, which was obtained by lyophilizing its 0.001M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with PTHrP.

**Immunogen:** PTHrP (15-34)-NH<sub>2</sub> (human), carrier free **Host:** Rabbit

### **Amino Acid Sequence of PTHrP (15-34)-NH<sub>2</sub> (human)<sup>1)</sup>**

AVSEHQLLHD KGKSIQDLRR RFFLHHLIAE IHTA-NH<sub>2</sub>

The amino acid sequences of N-terminal (15-34) of PTHrPs in mammals are 100% conserved.

**Product Form:** Lyophilized unpurified serum **Size:** 50 µL

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN<sub>3</sub> 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

**Suggested Working Dilution Range:** 1:1,000 (final dilution ~1:7,000) for radioimmunoassay; 1:200-1,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

**Specificity** (based on radioimmunoassay): PTHrP (15-34)-NH<sub>2</sub> 100%, PTHrP (1-34)-NH<sub>2</sub> 70%, PTHrP (1-19) 0%, PTHrP (7-34)-NH<sub>2</sub> 25%

**Positive Control** (immunohistochemistry): Rat mammary gland (lactation period) or mammary tumor tissue

**Species Tested:** Rat

### **REFERENCES:**

1) L.J. Suva, G.A. Winslow et al., A parathyroid hormone-related protein implicated in malignant hypercalcemia: cloning and expression. Science 237: 893-898, 1987

### **FOR RESEARCH LABORATORY USE ONLY**

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

#### **Distributor**



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

#### **Manufacturer**



2480-1 Awakura, Fujinomiya-shi,  
Shizuoka 418-0011, JAPAN  
Phone: +81-544-22-2771  
Fax: +81-544-22-2770  
URL : [www.yanaihara.co.jp](http://www.yanaihara.co.jp)  
E-mail : [ask@yanaihara.co.jp](mailto:ask@yanaihara.co.jp)