



**Anti Amyloid b Protein (1-40) (Human) Serum**  
**Cat. No. YII-YA010-EX      Lot No. 305100514**

**Description:** This antiserum was raised in a rabbit by immunization with a carrier free synthetic amyloid b protein (1-40) (human) peptide. The product vial contains 50 µL of the titled serum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with amyloid b protein (human).

**Immunogen:** Synthetic amyloid b protein (1-40) (human), carrier free **Host:** Rabbit

**Amino Acid Sequence of Amyloid b Protein (1-40) (human)<sup>1)</sup>:**

1 40

DAEFRHDSGY EVHHQKLVFF AEDVGSNKGA IIGLMVGGVV

**Product Form:** Lyophilized unpurified serum **Size:** 50 µL

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN<sub>3</sub> 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstituted, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

**Suggested Working Dilution Range:** 1:1,000 for immunohistochemistry. Optimal dilution should be determined by each laboratory for each application.

**Specificity** (based on enzyme immunoassay): Amyloid b protein (human) specific.

**REFERENCES:**

1) J. Kang, H.G. Lemaire et al., The precursor of Alzheimer's disease amyloid A4 protein resembles a cell-surface receptor. Nature 325:733-738, 1987

**FOR RESEARCH LABORATORY USE ONLY**

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

**Distributor**



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

**Manufacturer**



2480-1 Awakura, Fujinomiya-shi,  
Shizuoka 418-0011, JAPAN  
Phone: +81-544-22-2771  
Fax: +81-544-22-2770  
URL : [www.yanaihara.co.jp](http://www.yanaihara.co.jp)  
E-mail : [ask@yanaihara.co.jp](mailto:ask@yanaihara.co.jp)