

試験研究用となります。それ以外の用途のご使用はできません。

Catalog No. PEC-0003

TAT HIV-1 (48-60)

ペプチド名	TAT HIV-1 (48-60)
純度	≥95%
ペプチド配列	[H]GRKKRRQRRRPPQ[OH]
分子式	C70H131N35O16
分子量	1719.03 Da.
蛍光・標識等	—
Lot#	PP00000050-002
カウンターイオン	塩酸
推奨溶解溶媒	水
内容量	0.5mg (0.1mg×5)

使用・保管方法 -20℃以下で保管してください。長期保管する場合は-80℃以下を推奨します。
溶解したサンプルは短時間でご使用下さい。



コスモ・バイオ株式会社
COSMO BIO CO., LTD.

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル

URL : <http://http://www.cosmobio.co.jp/>

●営業部（お問合せ）

TEL : (03) 5632-9610

FAX : (03) 5632-9619

●札幌事業部（技術的なお問合せ）

TEL : (03) 5632-9744

FAX : (0134) 61-2295

For test research. It cannot be used for any other purpose.

Catalog No. PEC-0003

TAT HIV-1 (48-60)

Peptide name	TAT HIV-1 (48-60)
Purity	≥95%
Sequence	[H]GRKKRRQRRRPPQ[OH]
Molecular formula	C70H131N35O16
molecular weight	1719.03 Da.
fluorescence	–
Lot#	PP00000050-002
Counter ion	HCl
Dissolved solvent	Water
Quantity	0.5mg (0.1mg×5)

Storage Store below -20 °C. For longterm storage, recommend -80°C or below.
Please use the dissolved sample in a short time.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com