

For test research. It cannot be used for any other purpose.

Catalog No. RMP-0001

c-Myc Peptide

Peptide name	c-Myc Peptide
Purity	≥95%
Sequence	[H]EQKLISEEDL[OH]
Molecular formula	C ₅₁ H ₈₆ N ₁₂ O ₂₁
molecular weight	1,203.30 Da.
fluorescence	–
Lot#	PP00000045-001
Counter ion	TFA
Dissolved solvent	Water
Quantity	2.0 mg (1.0 mg×2)

Storage Store below -20 °C. For longterm storage, recommend -80°C or below.
Please use the dissolved sample in a short time.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
URL: www.cosmobiousa.com
Phone/FAX: (+1) 760-431-4600

試験研究用となります。それ以外の用途のご使用はできません。

Catalog No. RMP-0001

c-Myc Peptide

ペプチド名	c-Mycペプチド
純度	≥95%
ペプチド配列	[H]EQKLISEEDL[OH]
分子式	C ₅₁ H ₈₆ N ₁₂ O ₂₁
分子量	1,203.30 Da.
蛍光・標識等	—
Lot#	PP00000045-001
カウンターイオン	トリフルオロ酢酸
推奨溶解溶媒	水
内容量	2.0 mg (1.0 mg×2)

使用・保管方法 -20℃以下で保管してください。長期保管する場合は-80℃以下を推奨します。溶解したサンプルは短時間でご使用下さい。



コスモ・バイオ株式会社
COSMO BIO Co., LTD.

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル

URL:<http://www.cosmobio.co.jp/>

●営業部(お問合せ)

TEL:(03) 5632-9610

TEL:(03) 5632-9620

●札幌事業部(技術的なお問合せ)

FAX:(03) 5632-9619 TEL:(03) 5632-9744 FAX:(0134) 61-2295