



17(ALPHA), 20(BETA), 21-tri OH-Progesterone-HRP

Form Liquid
Volume 1 ml (1000TEST)

Application notes ELISA, RIA

Recommended use together **Antibody for ELISA:** Anti 17ALPHA, 20BETA, 21-tri OH-Progesterone Purified by Protein A (2000test for ELISA) 50µl/well use (Cat. no. FKA-340-E)
Antibody for RIA : Anti 17ALPHA, 20BETA, 21-tri OH-Progesterone (Cat. no.FKA-340)
Standard : 17A-20B-21TRI OH Progesterone (Cat.no.FKA-339ST)

Recommended dilutions ELISA : 50 µl/Well use
Optimal dilutions/concentrations should be determined by the end user.

Storage Maintain at 4°C-8°C in undiluted aliquots. Avoid freeze/thaw cycles

References The Cross reactivities of anti 17Alpha, 20Beta, 21-triOH-Pregnen-3-one-3-CMO-BSA(FKA-340E) by ELISA with 17Alpha, 20Beta, 21-triOH-Progesterone-3-MO-HRP (FKA-339)

O.D. (492/600nm) of B0 about 1.0
17Alpha, 20Beta, 21-triOH-Pregesterone dose of B0/2 about 2 ng/ml

Steroid

Cross Reaction

20Beta-S	100 %
Testosteron	1.10 %
Progesterone	1.00 %
Cortisol	0.90 %
17Alpha 20Alpha-diOH-P	0.80 %
20Alpha-S	0.28 %
Deoxycortisol(S)	0.14 %
Deoxycorticosterone(DOC)	0.07 %
17Alpha-OH-P	<0.01 %
17Alpha,20Beta-diOH-P	<0.01 %
Estradiol	<0.01 %



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619