Technical Support: tech-kobe@soyaku.co.jp


## Medicinal Chemistry Pharmaceutical Co., Ltd.

碠 Trans Genic Inc.

KZ622	mouse IgG1 isotype control antibody (Clone No.G1KD6)		
Туре	Monoclonal Antibody		Keyword
Immunogen	KLH		
Raised in	Mouse		
Myeloma	P3U1		
Clone number	G1KD6		
Isotype	lgG1, kappa		
Source	Serum-free medium		
Purification notes	ProteinG		Application
Cross Reactivity		WB	0.1-1.0 μg/mL
Label	Unlabeled	IHC	0.1-10 μg/mL
Formulation	Lyophilized powder	ICC	Not tested
Contents	100 μg	ELISA	Not tested
Reconstitution	Reconstitute at 100 μg/100 μL in DDW.	FCM	Not tested
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]	Neutralization	Not tested
Storage	This product is shipped at ambient temperature. Upon receipt, we recommend to store it immedeately below -20°C	IP	Not tested



KZ622



Sample: monkey hippocampus tissue lysate Primary antibody concentration: 1 µg/ml Secondary antibody: HRP-conjugated anti-mouse IgG

Note

Reference

## WARNING AND PRECAUTION

- 1. Not for diagnostic use. The safety and efficacy of product in diagnostic or other clinical uses has not been established.
- 2. Harmful by inhalation, in contact with skin and if swallowed. Do not breathe dust. Avoid contact with skin and eyes.
- 3. If contact with skin and eyes, wash all affected areas with large volume of water. If inhaled remove to fresh air. In severe case obtain medical attention.
- 4. Wash hand thoroughly after handling the product.5. Do not use this product if container is broken or some contaminants are detected.
- 6. When preserving the product, Close the container, ensure it does not fall aside or
- 7. Dispose of the container and expired reagents in accordance with federal, state and local government regulations.
- 8. Do not use the container and accessories of the product for other purpose.

## 取り扱い上の注意

- この添付文書をよく読んでから使用して下さい。
- 1. 本品は、研究用試薬であり、医薬品その他の目的にはご使用になれません。
- 2. 取り扱い中は皮膚、粘膜、着衣に触れたり、目に入らないように適切な措置を行って下さ
- 3. 誤って目や口に入った場合は、水で十分に洗い流すなどの応急処置を行い、必要があ れば医師の手当を受けて下さい。
- 4.取り扱い後には手洗いを十分に行って下さい。
- 5. 容器の破損、異物混入等異常が認められた物は使用しないで下さい。
- 6. 試薬を保管する場合は、蓋をし、転倒落下防止を確実にし、指定の貯蔵方法で保管して
- 7. 使用後の容器は、廃棄物に関する規定に従って処理して下さい。
- 8. 容器、付属品等の他目的への転用は保証できません。