



BPO-BGG

Form	Lyophilized protein powder Benzylpenicilloyl Bovine γ -Globulin
Volume	Protein: 20mg (As average of molecular weight: 170,000)
Storage	Store below 4°C (below -20°C for prolonged storage or after dissolved)
Concentration	Add distilled water. It will become 5mM boric acid buffered saline [BBS] , 0.9% NaCl solution.

Note This product is lyophilized product containing totally 20mg proteins per vial. Therefore, conjugation ratio BPO groups/molecule, and volume of distilled water to reconstitute is different by production lots.

For detailed information, visit the website below

https://www.cosmobio.com/products/lsl_lg0577_1.html

Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG

LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-2334	Mite-Df Feces AG
LG-2444	Mite-Dp Feces AG
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

Adjuvant

LG-6000	ALUM
---------	------

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: info@cosmobiousa.com

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

BPO-BGG

I. 内容

lot No. 812120

性 状	Benzylpenicilloyl Bovine γ -Globulin 凍結乾燥粉末
規 格	蛋白質量: 20mg BPO基数: 1分子あたり 5.3個 (平均分子量 170,000 としての計算値) 1.0mlの精製水を加えることで5mM硼酸緩衝溶液 (pH 8.0) 0.9% NaCl 溶液になります。

II. 保存上の注意

4℃以下 (長期間又は溶液の場合は-20℃以下) で保存して下さい。

製造元

総発売元

LSL 株式会社 エル・エス・エル

コスモ・バイオ株式会社

新製品案内

- ☆Mite Extract
カタログNo. LG-5119 (LSL)
- ☆Cedar Pollen Extract
カタログNo. LG-5229 (LSL)
- ☆BPO-BSA
カタログNo. LG-0317 (LSL)
- ☆BPO-BGG
カタログNo. LG-0577 (LSL)
- ☆DNP-BSA (30)
カタログNo. LG-3017 (LSL)
- ☆DNP-KLH
カタログNo. LG-0089 (LSL)