



Mite Extract-Tp

Form	Lyophilized protein powder Mite-Tp Crude Extract This product extracted from House Dust Mites, <u>Tyrophagus putrescentiae</u> (Schrank) *. * Mites were homogenized with 10-fold 50 mM Phosphate buffer (pH7.2) and the homogenate were stirred at 4°C for 24 hours. And centrifugal supernatant dialyzed by PBS.
Volume	Protein: 10mg
Storage	Store below 4°C (below -20°C for prolonged storage or after dissolved)
Concentration	Add distilled water. It will become 5mM phosphate buffered saline [PBS] , 0.9% NaCl solution.
Note	This product is lyophilized product containing totally 10mg proteins per vial. Therefore, conjugation ratio DNP groups/molecule, and volume of distilled water to reconstitute is different by production lots. For detailed information, visit the website below https://www.cosmobio.com/products/lsl_lg5559_1.html

Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG

LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-2334	Mite-Df Feces AG
LG-2444	Mite-Dp Feces AG
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

Adjuvant

LG-6000	ALUM
---------	------

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Mite Extract-Tp

I. 内容

lot No. 784041

性 状	凍結乾燥粉末 (Mite-Tp粗抽出物) House dust mites, <u>Tyrophagus putrescentiae</u> (Schrank) ケナガコナダニ (コナダニ科) の成虫から抽出*
規 格	蛋白質量: 10mg 3.3mlの精製水を加える事により5mM硼酸緩衝溶液 (pH 8.0) 0.9% NaCl 溶液になります。

* 虫体を10倍量の50mMリン酸緩衝液 (pH 7.2) でホモジェナイズ
、24時間攪拌 (4℃)、遠心上澄を濃縮後BBSで透析。

II. 保存上の注意

4℃以下 (長期間又は溶液の場合は-20℃以下) で保存して下さい。

製造元

総発売元

LSL 株式会社 エル・エス・エル

コスモ・バイオ株式会社

新製品案内

☆Mite Extract	カタログNo. LG-5119 (LSL)
☆Cedar Pollen Extract	カタログNo. LG-5229 (LSL)
☆Mite Extract-Df	カタログNo. LG-5339 (LSL)
☆Mite Extract-Dp	カタログNo. LG-5449 (LSL)
☆Mite Extract-Tp	カタログNo. LG-5559 (LSL)
☆Anti DNP (Mono)	カタログNo. LM-0909 (LSL)
☆Anti Mite Extract	カタログNo. LB-5199 (LSL)