



Hinoki Cypress, Pollen Crude Extract

Form	Lyophilized protein powder Hinoki Pollen Crude Extract
	This product extracted from Hinoki Cypress, <u>Chamaecyparis obtusa</u> Pollen*. *Pollens were stirred with 20 volumes of 0.125M NaHCO ₃ (pH8) at 4°C for 16 hours, homogenized, centrifuged, and the supernatant dialyzed against BBS.
Volume	Protein: 2mg
Concentration	Add 2.0mL of distilled water. It will become 5mM Borate buffered saline [BBS] (pH 8.0), 0.9% NaCl solution.
Storage	Store below 4°C (below -20°C for prolonged storage or after dissolved)

Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG
LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

Adjuvant

LG-6000	ALUM
---------	------

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Hinoki Pollen Extract

I. 内容

lot No. 737092

性 状	凍結乾燥粉末 (Hinoki Pollen粗抽出物) Hinoki Cypress, <u>Chamaecyparis obtusa</u> の花粉から抽出*
規 格	蛋白質量: 2mg 2.0mlの精製水を加える事により5mM硼酸緩衝溶液 (pH8.0) 0.9% NaCl [BBS] 溶液になります。

*花粉を20倍量の 0.125M NaHCO₃ (pH8) で16時間攪拌 (4℃) 後
Homogenize、遠心上澄を濃縮後BBSで透析 (Visking tube, 4℃)。

II. 保存上の注意

4℃以下 (長期間又は溶液の場合は-20℃以下) で保存して下さい。

製造元

総発売元

LSL 株式会社 エル・エス・エル

コスモ・バイオ株式会社

新製品案内

☆Mite Extract-Df	カタログNo. LG-5339 (LSL)
☆Mite Extract-Dp	カタログNo. LG-5449 (LSL)
☆Mite Extract-Tp	カタログNo. LG-5559 (LSL)
☆Hinoki Pollen Extract	カタログNo. LG-5779 (LSL)
☆Ragweed Pollen Extract	カタログNo. LG-5779 (LSL)
☆百日咳死菌 (FTPP) 懸濁液	カタログNo. LG-7000 (LSL)