



Ragweed Pollen Extract

Form
Lyophilized protein powder
Ragweed Pollen Crude Extract

This product produced from Short Ragweed, Ambrosia artemisiifolia*.
*Acetone treated pollen was stirred with 20 volumes of 0.125M sodium bicarbonate (pH 8) at 4°C for 16 hours, centrifuged, and the supernatant dialyzed against PBS.

Volume
Protein: 5mg

Concentration
Add 1.0mL of distilled water. It will become 5mM Borate buffered saline [BBS] (pH 8.0), 0.9% NaCl solution.

Storage
Store below 4°C
(below -20°C for prolonged storage or after dissolved)

Note
This product is lyophilized product containing totally 5mg proteins per vial. Therefore, volume of distilled water to reconstitute is different by production lots. For detailed information, visit the website below
https://www.cosmobio.com/products/lsl_lg5889_1.html

Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG

LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-2334	Mite-Df Feces AG
LG-2444	Mite-Dp Feces AG
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

Adjuvant

LG-6000	ALUM
---------	------

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Ragweed Pollen Extract

I. 内容

lot No. 726091

性 状	凍結乾燥粉末 (Ragweed Pollen粗抽出物) Short Ragweed, <u>Ambrosia artemisiifolia</u> の花粉から抽出 *
規 格	蛋白質量: 5mg 1.0mlの精製水を加える事により5mM硼酸緩衝溶液 (pH8.0) 0.9% NaCl溶液になります。

* アセトン処理した花粉を20倍量の 0.125M NaHCO₃ (pH8) で16時間攪拌
(4℃)、遠心上澄を濃縮後BBSで透析。

II. 保存上の注意

4℃以下 (長期間又は溶液の場合は-20℃以下) で保存して下さい。

製造元

総発売元

LSL 株式会社 エル・エス・エル

コスモ・バイオ株式会社

新製品案内

☆Mite Extract
カタログNo. LG-5119 (LSL)
☆Mite Extract-Df
カタログNo. LG-5339 (LSL)
☆Mite Extract-Dp
カタログNo. LG-5449 (LSL)
☆Mite Extract-Tp
カタログNo. LG-5559 (LSL)
☆Cedar Pollen Extract
カタログNo. LG-5229 (LSL)
★Cedar Pollen Extract-Cj
カタログNo. LG-5280 (LSL)
☆Anti Mite Extract
カタログNo. LB-5199 (LSL)