



Mite Extract-Dp

Form Lyophilized protein powder
Mite-Dp Crude Extract

This product extracted from House Dust Mites, *Dermatophagoides pteronyssinus* (Trouessari)*.

*Mites were homogenized with 10-fold 50 mM Phosphate buffer (pH7.2) and the homogenate were stirred at 4°C for 24 hours. And centrifugal supernatant dialyzed 4 times by PBS (Visking tube, 4°C).

Volume Protein: 10mg

Storage Store below 4°C
(below -20°C for prolonged storage or after dissolved)

Concentration Add distilled water.
It will become phosphate buffered saline [PBS] (Dulbeco), solution.

Note This product is lyophilized product containing totally 10mg proteins per vial. Therefore, volume of distilled water to reconstitute is different by production lots. For detailed information, visit the website below
https://www.cosmobio.com/products/lsl_lg5449_1.html

Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG

LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-2334	Mite-Df Feces AG
LG-2444	Mite-Dp Feces AG
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

Adjuvant

LG-6000	ALUM
---------	------

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Mite Extract-Dp

I. 内容

lot No. 136101

性 状	凍結乾燥粉末 (Mite-Dp粗抽出物) House dust mites, <i>Dermatophagoides pteronyssinus</i> (Trouessari) ヤケヒョウヒダニ (チリダニ科) の成虫から抽出*
規 格	蛋白質量: 10mg 3.6mlの精製水を加える事によりPBS (Dulbecco, pH7.4) 溶液になります。

* 虫体を10倍量の50mMリン酸緩衝液 (pH7.2) でホモジェナイズ、24時間攪拌 (4℃)、遠心上澄を濃縮後PBSで透析 (Visking tube, 4℃)。既知アレルゲンの含有量はGroup1が約90μg/vial、Group2が約320μg/vial 程度です。

II. 保存上の注意

4℃以下 (長期間又は溶液の場合は-20℃以下) で保存して下さい。

製造元

総発売元

LSI株式会社 エル・エス・エル



コスモ・バイオ株式会社

※LSL社の抽出物は、蛋白質抽出時の濃度がロットによって異なります。ロット毎に加える精製水の量も異なりますので、商品に添付されているデータシートをご確認ください。

新製品案内

LG-2334	Mite-Df Feces AG (ダニ-Df由来排泄物抗原)
LG-2444	Mite-Dp Feces AG (ダニ-Dp由来排泄物抗原)
LG-5999	Wormwood Pollen Extract
LG-7000	FTPP Suspension (百日咳死菌懸濁液)
LG-8000	FTMB Suspension (結核死菌懸濁液)
LB-7111	Anti (MAP) Mite Der fl
LB-7103	Anti (MAP) Mite Group 1
LB-7204	Anti (MAP) Mite Group 2