

## Aggrecan Antibody to N-terminal neoepitope FFGV Mouse monoclonal Antibody

Catalog Number 1042004

<b>INTRODUCTION</b>	This fragment is rapidly released from the tissue when MMP catabolism of aggrecan occurs and has been identified in synovial fluid samples from patients with degenerative joint diseases.
<b>IMMUNOGEN</b>	Synthetic peptide: FFGVGGE.
<b>CLONE</b>	BC-14
<b>HOST</b>	Mouse
<b>MYELOMA</b>	x63-Ag8.653
<b>ISOTYPE</b>	IgG2a
<b>LIGHT CHAIN TYPE</b>	kappa
<b>SPECIFICITY</b>	Recognizes the N-terminal neoepitope sequence (FFGV....) generated at the “MMP cleavage site” after MMP catabolism in the interglobular domain of aggrecan between amino acids ..PEN341 and 342FFG.. (Human sequence enumeration).  This antibody cross-reacts with Human, Rat, Cow, Guinea pig, Horse, Pig.
<b>PURITY</b>	Affinity purified on protein G

*For research use only.*

North America: 1-888-USMDBIO | International: +41-44 986 2628 | [products@mdbiosciences.com](mailto:products@mdbiosciences.com) | [www.mdbioproductions.com](http://www.mdbioproductions.com)



IN30140.00

<b>STORAGE BUFFER</b>	PBS, no preservatives
<b>FORM</b>	Liquid, 1 mL/vial
<b>CONCENTRATION</b>	0.1 mg/mL
<b>APPLICATION</b>	<ul style="list-style-type: none"><li>• Western-Blotting Suggested dilution: 1:100 Detects a a band of approximately 60 - 150 kDa</li><li>• ELISA</li><li>• IHC</li></ul>
<b>TECHNICAL NOTES</b>	<p>This antibody should work in IHC on formalin- or paraformaldehyde-fixed paraffin embedded sections as well as either alcohol-fixed frozen sections or un-fixed snap-frozen sections.</p> <p>Samples are usually deglycosylated using 0.01 Units Chondroitinase ABC (Sigma), 0.01 Units Keratanase (Seikagaku) and 0.0001 Units Keratanase II (Seikagaku) per 10µg S-GAG of non-deglycosylated aggrecan for optimal epitope recognition in SDS-PAGE and immunohistochemistry. (1, 2).</p>
<b>STORAGE</b>	The antibody is stable until the expiry date given on the label if stored at -20 °C. Repeated freezing and thawing should be avoided. Aliquoting is recommended.
<b>EXPIRATION</b>	See Vial Label
<b>REFERENCES</b>	<ol style="list-style-type: none"><li>1. Little CB et al. Matrix metalloproteinases are involved in C-terminal and interglobular domain processing of cartilage aggrecan in late stage cartilage degradation. Matrix Biol 21:271-88 (2002).</li><li>2. Caterson B et al. Mechanisms involved in cartilage proteoglycan catabolism. Matrix Biol 19:333-44 (2000)</li></ol>

*For research use only.*

North America: 1-888-USMDBIO | International: +41-44 986 2628 | products@mdbiosciences.com | www.mdbioproducts.com

