

STABLE LIQUID ALKALINE PHOSPHATASE
INSOLUBLE SUBSTRATES

FAST RED

ASLB-100 100 mls
ASLB-500 500 mls
ASLB-1000 1000 mls

Two components.

Stable for 18 months at 2-8° C.

Used for tissue sections or immunoblotting. Uses LB salt which produces a brilliant red color at reaction site. Easy to use. Simply mix the 2 liquid components 1:1 and apply.

Moss, Inc., P.O. Box 189, Pasadena, Maryland 21123-0189, U.S.A.

Headquarters Phone: (800) 932-6677, (410) 768-3442 Fax: (410) 768-3971

Southern California office: Phone: (949) 552-1652 Fax: (949) 552-6807

Northern California office: Phone: (650) 343-4595 Fax: (650) 343-8420

