

# JORDAN COMBS

COMBS

## JSB-30 - COMBS

| CATALOG# | DESCRIPTION      | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|------------------|------------|--------------------------|
| CSM-0808 | 0.8mm x 8 tooth  | 5.8mm      | 4.7µl                    |
| CSM-0813 | 0.8mm x 13 tooth | 4.0mm      | 3.2µl                    |
| CSM-0818 | 0.8mm x 18 tooth | 2.3mm      | 1.9µl                    |
| CSM-1508 | 1.5mm x 8 tooth  | 5.8mm      | 8.7µl                    |
| CSM-1513 | 1.5mm x 13 tooth | 4.0mm      | 6.0µl                    |
| CSM-1518 | 1.5mm x 18 tooth | 2.3mm      | 3.5µl                    |

## JSB-302 - COMBS

| CATALOG# | DESCRIPTION      | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|------------------|------------|--------------------------|
| CSD-0816 | 0.8mm x 16 tooth | 5.8mm      | 4.7µl                    |
| CSD-0827 | 0.8mm x 27 tooth | 4.0mm      | 3.2µl                    |
| CSD-0836 | 0.8mm x 36 tooth | 2.3mm      | 1.9µl                    |
| CSD-1516 | 1.5mm x 16 tooth | 5.8mm      | 8.7µl                    |
| CSD-1522 | 1.5mm x 22 tooth | 5.0mm      | 7.5µl                    |
| CSD-1527 | 1.5mm x 27 tooth | 4.0mm      | 6.0µl                    |
| CSD-1533 | 1.5mm x 33 tooth | 3.0mm      | 4.5µl                    |
| CSD-1536 | 1.5mm x 36 tooth | 2.3mm      | 3.5µl                    |

## JSB-96 - COMBS

| CATALOG#  | DESCRIPTION                | WELL WIDTH   | SAMPLE VOLUME PER MM GEL |
|-----------|----------------------------|--------------|--------------------------|
| CB96-0825 | 0.8mm x 25 tooth, MCP**    | 5.8mm        | 4.7µl                    |
| CB96-0840 | 0.8mm x 40 tooth           | 4.0mm        | 3.2µl                    |
| CB96-0850 | 0.8mm x 50 tooth, MCP**    | 2.9mm        | 2.3µl                    |
| CB96-1525 | 1.5mm x 25 tooth, MCP**    | 5.8mm        | 8.7µl                    |
| CB96-1540 | 1.5mm x 40 tooth           | 4.0mm        | 6.0µl                    |
| CB96-1550 | 1.5mm x 50 tooth, MCP**    | 2.9mm        | 4.4µl                    |
| CB96-3025 | 3.0mm x 25 tooth, MCP**    | 5.8mm        | 17.4µl                   |
| CB96-3040 | 3.0mm x 40 tooth           | 4.0mm        | 12.0µl                   |
| CB96-15P  | 1.5mm, 2 markers, 1 sample | 6.4, 204.0mm | 9.6, 306.0µl             |
| CB96-30P  | 3.0mm, 2 markers, 1 sample | 6.4, 204.0mm | 19.2, 612.0µl            |

## JSB-75 - COMBS

| CATALOG#  | DESCRIPTION      | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|-----------|------------------|------------|--------------------------|
| CB75-0808 | 0.8mm x 8 tooth  | 9.9mm      | 7.9µl                    |
| CB75-0813 | 0.8mm x 13 tooth | 5.8mm      | 4.7µl                    |
| CB75-0822 | 0.8mm x 22 tooth | 4.0mm      | 16.0µl                   |
| CB75-0830 | 0.8mm x 30 tooth | 2.5mm      | 2.0µl                    |
| CB75-1508 | 1.5mm x 8 tooth  | 9.9mm      | 14.9µl                   |
| CB75-1513 | 1.5mm x 13 tooth | 5.8mm      | 8.7µl                    |
| CB75-1522 | 1.5mm x 22 tooth | 4.0mm      | 6.0µl                    |
| CB75-1530 | 1.5mm x 30 tooth | 2.5mm      | 3.7µl                    |
| CB75-3008 | 3.0mm x 8 tooth  | 9.9mm      | 29.72µl                  |
| CB75-3013 | 3.0mm x 13 tooth | 5.8mm      | 17.5µl                   |
| CB75-3022 | 3.0mm x 22 tooth | 4.0mm      | 12.0µl                   |

## JSB-120 - COMBS

| CATALOG#   | DESCRIPTION      | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|------------|------------------|------------|--------------------------|
| CB120-0812 | 0.8mm x 12 tooth | 10.0mm     | 8.1µl                    |
| CB120-0821 | 0.8mm x 21 tooth | 5.8mm      | 4.7µl                    |
| CB120-0836 | 0.8mm x 36 tooth | 4.0mm      | 3.2µl                    |
| CB120-0848 | 0.8mm x 48 tooth | 2.5mm      | 2.0µl                    |
| CB120-1512 | 1.5mm x 12 tooth | 10.0mm     | 15.2µl                   |
| CB120-1521 | 1.5mm x 21 tooth | 5.8mm      | 8.7µl                    |
| CB120-1536 | 1.5mm x 36 tooth | 4.0mm      | 6.0µl                    |
| CB120-1548 | 1.5mm x 48 tooth | 2.5mm      | 3.7µl                    |
| CB120-3012 | 3.0mm x 12 tooth | 10.0mm     | 30.5µl                   |
| CB120-3021 | 3.0mm x 21 tooth | 5.8mm      | 17.5µl                   |
| CB120-3036 | 3.0mm x 36 tooth | 4.0mm      | 12.0µl                   |

MCP\*\* - Multichannel Pipette Compatible

## JVD-80 - COMBS

| CATALOG# | DESCRIPTION                | WELL WIDTH   | VOLUME PER MM WELL DEPTH |
|----------|----------------------------|--------------|--------------------------|
| CP-0406  | 0.4mm x 6 tooth            | 25.4mm       | 10.1µl                   |
| CP-0412  | 0.4mm x 12 tooth           | 9.9mm        | 4.0µl                    |
| CP-0418  | 0.4mm x 18 tooth           | 5.8mm        | 2.3µl                    |
| CP-0427  | 0.4mm x 27 tooth           | 4.0mm        | 1.6µl                    |
| CP-0806  | 0.8mm x 6 tooth            | 25.4mm       | 20.3µl                   |
| CP-0812  | 0.8mm x 12 tooth           | 9.9mm        | 7.9µl                    |
| CP-0818  | 0.8mm x 18 tooth           | 5.8mm        | 4.6µl                    |
| CP-0827  | 0.8mm x 27 tooth           | 4.0mm        | 3.2µl                    |
| CP-1506  | 1.5mm x 6 tooth            | 25.4mm       | 38.1µl                   |
| CP-1512  | 1.5mm x 12 tooth           | 9.9mm        | 14.9µl                   |
| CP-1518  | 1.5mm x 18 tooth           | 5.8mm        | 8.7µl                    |
| CP-1527  | 1.5mm x 27 tooth           | 4.0mm        | 6.0µl                    |
| CP-3006  | 3.0mm x 6 tooth            | 25.4mm       | 76.2µl                   |
| CP-3012  | 3.0mm x 12 tooth           | 9.9mm        | 29.7µl                   |
| CP-3018  | 3.0mm x 18 tooth           | 5.8mm        | 17.4µl                   |
| CP-3027  | 3.0mm x 27 tooth,          | 4.0mm        | 12.0µl                   |
| CP-04P   | 0.4mm, 2 markers, 1 sample | 6.3, 152.4mm | 2.5, 61.0µl              |
| CP-08P   | 0.8mm, 2 markers, 1 sample | 6.3, 152.4mm | 5.0, 121.9µl             |
| CP-15P   | 1.5mm, 2 markers, 1 sample | 6.3, 152.4mm | 9.5, 228.6µl             |
| CP-30P   | 3.0mm, 2 markers, 1 sample | 6.3, 152.4mm | 18.9, 457.2µl            |
| PF-26P   | 0.4mm x 26 tooth*          | 3.0mm        | 1.2µl                    |
| PF-36P   | 0.4mm x 32 tooth*          | 2.5mm        | 1.0µl                    |
| PF-46P   | 0.4mm x 46 tooth*          | 2.0mm        | 0.8µl                    |

## JQT-182 - COMBS

| CATALOG# | DESCRIPTION             | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|-------------------------|------------|--------------------------|
| CQ-0418  | 0.4mm x 18 tooth        | 13.6mm     | 5.44µl                   |
| CQ-0432  | 0.4mm x 32 tooth        | 5.8mm      | 2.32µl                   |
| CQ-0454  | 0.4mm x 54 tooth        | 4.0mm      | 1.6µl                    |
| PF-64Q   | 0.4mm x 64 picket fence | 1.9mm      | 0.8µl                    |
| PF-80Q   | 0.4mm x 80 picket fence | 1.9mm      | 0.8µl                    |



IBI ANALYTICAL COMBS



IBI PREPARATORY COMBS