

IB50000 - QSH LAB-PAL COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|----------------------------|-------------|--------------------------|
| IB50040 | 1.0mm x 5 tooth | 6.0mm | 6.0μl |
| IB50050 | 1.0mm x 10 tooth | 2.6mm | 2.6μl |
| IB50080 | 1.5mm x 5 tooth | 6.0mm | 9.0μl |
| IB50090 | 1.5mm x 10 tooth | 2.6mm | 3.9μl |
| IB50060 | 1.0mm, 0 marker, 1 sample | 44.0mm | 44.0μl |
| IB50070 | 1.0mm, 1 marker, 1 sample | 6.0, 34.6mm | 6.0, 34.6μl |
| IB50100 | 1.5mm, 0 marker, 1 sample | 44.0mm | 66.0μl |
| IB50110 | 1.5mm, 1 marker, 1 samples | 6.0, 34.6mm | 9.0, 52.0μl |

IB51000 - QS-710 COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|---------------------------|-------------|--------------------------|
| IB51040 | 1.0mm x 8 tooth | 4.7mm | 4.7μl |
| IB51045 | 1.0mm x 12 tooth | 2.8mm | 2.8μl |
| IB51050 | 1.0mm x 15 tooth | 2.2mm | 2.2μl |
| IB51055 | 1.5mm x 8 tooth | 4.7mm | 7.0μl |
| IB51060 | 1.5mm x 12 tooth | 2.8mm | 4.0μl |
| IB51065 | 1.5mm x 15 tooth | 2.2mm | 3.3μl |
| IB51100 | 1.5mm x 3 tooth sets | 3.0, 9.5mm | 4.5, 14.25μl |
| IB51101 | 1.5mm x 2 tooth sets | 3.0, 9.5mm | 4.5, 14.25μl |
| IB51046 | 2.0mm x 12 tooth | 2.8mm | 6.0μl |
| IB51090 | 2.0mm x 10 tooth | 4.0mm | 8.0μl |
| IB51095 | 3.0mm x 6 tooth | 5.8mm | 17.3μl |
| IB51070 | 1.0mm, 0 marker, 1 sample | 62.6mm | 62.6μl |
| IB51075 | 1.0mm, 1 marker, 1 sample | 4.7, 54.7mm | 4.7, 54.7μl |
| IB51080 | 1.5mm, 0 marker, 1 sample | 62.6mm | 94.0μl |
| IB51085 | 1.5mm, 1 marker, 1 sample | 4.7, 54.7mm | 7.0, 82.0μl |

IB53000 / IB52000 - MP-1015 / MPH COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|-----------------------------|-------------|--------------------------|
| IB52060 | 1.0mm x 10 tooth | 6.3mm | 6.3μl |
| IB52090 | 1.0mm x 16 tooth | 3.5mm | 3.5μl |
| IB52160 | 1.0mm x 20 tooth | 2.6mm | 2.6μl |
| IB52165 | 1.5mm x 20 tooth | 2.6mm | 4.5μl |
| IB52155 | 1.5mm x 22 tooth | 2.6mm | 3.9μl |
| IB53100 | 1.5mm x 5 tooth sets | 3.0, 9.5mm | 4.5, 14.25μl |
| IB53101 | 1.5mm x 4 tooth sets | 3.0, 9.5mm | 4.5, 14.25μl |
| IB52070 | 2.0mm x 10 tooth | 6.3mm | 12.6μl |
| 8920316 | 2.0mm x 11 tooth, MCP** | 6.3mm | 12.6μl |
| IB53070 | 2.0mm x 13 tooth | 4.0mm | 8.0μl |
| IB52100 | 2.0mm x 16 tooth | 3.5mm | 7.0μl |
| IB52170 | 2.0mm x 20 tooth | 2.6mm | 5.2μl |
| IB52080 | 3.0mm x 10 tooth | 6.3mm | 18.9μl |
| 8695108 | 3.0mm x 11 tooth, MCP** | 6.3mm | 18.9μl |
| IB52110 | 3.0mm x 16 tooth | 3.5mm | 10.5μl |
| IB52180 | 3.0mm x 20 tooth | 2.6mm | 7.8μl |
| IB52120 | 1.0mm, 2 markers, 1 sample | 2.9, 81.8mm | 2.9, 82.0μl |
| IB52130 | 2.0mm, 2 markers, 1 sample | 2.9, 81.8mm | 5.8, 163.0μl |
| IB53080 | 3.0mm, 2 markers, 1 sample | 2.9, 81.8mm | 8.7, 245.0μl |
| IB52140 | 1.0mm, 2 markers, 2 samples | 2.9, 39.3mm | 2.9, 39.0μl |
| IB52150 | 2.0mm, 2 markers, 2 samples | 2.9, 39.3mm | 5.8, 78.0μl |
| IB53090 | 3.0mm, 2 markers, 2 samples | 2.9, 39.3mm | 8.7, 117.0μl |

IB56000 / IB54000 - HR-2025 / HRH COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|-------------------------|------------|--------------------------|
| IB54180 | 1.0mm x 10 tooth | 16.6mm | 16.6μl |
| IB56050 | 1.0mm x 12 tooth | 13.1mm | 13.1μl |
| IB56055 | 1.0mm x 15 tooth | 9.8mm | 9.8μl |
| IB54060 | 1.0mm x 20 tooth | 6.5mm | 6.5μl |
| IB54090 | 1.0mm x 30 tooth | 4.0mm | 4.0μl |
| IB56100 | 1.0mm x 40 tooth | 3.2mm | 3.2μl |
| IB56101 | 1.0mm x 44 tooth, MCP** | 2.9mm | 2.9μl |
| IB54190 | 2.0mm x 10 tooth | 16.6mm | 49.8μl |
| IB56115 | 2.0mm x 12 tooth | 13.1mm | 26.2μl |

IB56000 / IB54000 - HR-2025 / HRH COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|-----------------------------|--------------|--------------------------|
| IB56125 | 2.0mm x 15 tooth | 9.8mm | 19.6μl |
| IB54070 | 2.0mm x 20 tooth | 6.5mm | 13.0μl |
| 8091100 | 2.0mm x 22 tooth, MCP** | 6.3mm | 12.6μl |
| IB54100 | 2.0mm x 30 tooth | 4.0mm | 8.0μl |
| IB56135 | 2.0mm x 40 tooth | 3.2mm | 6.4μl |
| IB56102 | 2.0mm x 44 tooth, MCP** | 2.9mm | 5.8μl |
| IB54080 | 3.0mm x 20 tooth | 6.5mm | 19.5μl |
| 8093767 | 3.0mm x 22 tooth, MCP** | 6.3mm | 18.9μl |
| IB54110 | 3.0mm x 30 tooth | 4.0mm | 12.0μl |
| IB56140 | 3.0mm x 40 tooth | 3.2mm | 9.6μl |
| IB56103 | 3.0mm x 44 tooth, MCP** | 2.9mm | 8.7μl |
| IB54120 | 1.0mm, 2 markers, 1 sample | 4.0, 184.0mm | 4.0, 184.0μl |
| IB54140 | 1.0mm, 2 markers, 2 samples | 4.0, 91.0mm | 4.0, 91.0μl |
| IB54130 | 2.0mm, 2 markers, 1 sample | 4.0, 184.0mm | 8.0, 368.0μl |
| IB54150 | 2.0mm, 2 markers, 2 samples | 4.0, 91.0mm | 8.0, 182.0μl |
| IB56145 | 3.0mm, 2 markers, 1 sample | 4.0, 184.0mm | 12.0, 552.0μl |
| IB56150 | 3.0mm, 2 markers, 2 samples | 4.0, 91.0mm | 12.0, 273.0μl |

IB57000 - HR2525 COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | SAMPLE VOLUME PER MM GEL |
|----------|-------------------------|------------|--------------------------|
| IB57125 | 1.0mm x 25 tooth, MCP** | 5.8mm | 5.8μl |
| IB57150 | 1.0mm x 50 tooth, MCP** | 2.9mm | 2.9μl |
| IB57225 | 2.0mm x 25 tooth, MCP** | 5.8mm | 11.6μl |
| IB57250 | 2.0mm x 50 tooth, MCP** | 2.9mm | 5.8μl |
| IB57325 | 3.0mm x 25 tooth, MCP** | 5.8mm | 17.5μl |
| IB57350 | 3.0mm x 50 tooth, MCP** | 2.9mm | 8.7μl |
| IB57425 | 1.5mm x 25 tooth, MCP** | 5.8mm | 8.7μl |
| IB57450 | 1.5mm x 50 tooth, MCP** | 2.9mm | 4.4μl |

IB80000 - STS-45i COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | VOLUME PER MM WELL DEPTH |
|----------|--------------------------|------------|--------------------------|
| IB80040 | 0.2mm x 32 tooth | 8.0mm | 1.6μl |
| IB80050 | 0.2mm x 64 tooth | 3.0mm | 0.6μl |
| IB80100 | 0.4mm x 32 tooth | 8.0mm | 3.2μl |
| IB80110 | 0.4mm x 64 tooth | 3.0mm | 1.2μl |
| IB80250 | 0.2mm x 64 sharktooth | 3.0mm | 0.6μl |
| IB80355 | 0.2mm x 64&96 sharktooth | 3.0, 2.0mm | 0.2, 0.4μl |
| IB80310 | 0.4mm x 64 sharktooth | 3.0mm | 1.2μl |
| IB80350 | 0.2mm x 96 sharktooth | 2.0mm | 0.4μl |
| IB80380 | 0.4mm x 96 sharktooth | 2.0mm | 0.8μl |
| IB80385 | 0.4mm x 64&96 sharktooth | 3.0, 2.0mm | 1.2, 0.8μl |

IB62000 - VCV COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | VOLUME PER MM WELL DEPTH |
|----------|----------------------------|--------------|--------------------------|
| IB62045 | 0.8mm x 12 tooth | 10.0mm | 8μl |
| IB62105 | 0.8mm x 20 tooth | 5.0mm | 4μl |
| IB62055 | 1.5mm x 12 tooth | 10.0mm | 15μl |
| IB62115 | 1.5mm x 20 tooth | 5.0mm | 7.5μl |
| IB62065 | 3.0mm x 12 tooth | 10.0mm | 30μl |
| IB62125 | 3.0mm x 20 tooth | 5.6mm | 15μl |
| IB62345 | 0.8mm, 2 markers, 2 sample | 6.3, 72.0mm | 5.0, 57.6μl |
| IB62355 | 1.5mm, 2 markers, 2 sample | 6.3, 72.0mm | 9.5, 108.0μl |
| IB62365 | 3.0mm, 2 markers, 1 sample | 6.3, 147.0mm | 18.9, 441.0μl |
| IB62375 | 3.0mm, 2 markers, 2 sample | 6.3, 72.0mm | 18.9, 216.0μl |

IB94000 - MAGELIN COMBS

| CATALOG# | DESCRIPTION | WELL WIDTH | VOLUME PER MM WELL DEPTH |
|----------|----------------------------|------------|--------------------------|
| IB92030 | 0.8mm x 10 tooth | 5.9mm | 28.1μl |
| IB92031 | 0.8mm x 12 tooth | 4.2mm | 20.3μl |
| IB92032 | 0.8mm x 15 tooth | 3.4mm | 16.2μl |
| IB92033 | 0.8mm, 0 markers, 1 sample | 7.0mm | 336.0μl |