



COSMO BIO CO., LTD.

Inspiration for Life Science

Agrisera 
Antibodies for plant sciences

total pre-immune IgY
AS03 027

Product information

Product number: AS03 027
Lot Number: 0403
Raised in: Hen
Purity: Total IgY
Quantity: 100 mg
Concentration: 22 mg/ml
in PBS pH 8
Antibody form: Liquid

Storage instructions: 4°C
*Preparation contains sodium azide
(0.02 %) as inhibitor of bacterial
growth.*

DO NOT FREEZE!

Background

Total pre-immune IgY is purified from egg yolks of hens not immunized with any antigens.

Please, remember to spin tubes briefly prior to opening them to avoid any losses that might occur from liquid or lyophilized material adhering to the cap or sides of the tubes.

Antibodies are intended for the research use only not for diagnostic or therapeutic use.

Distributor



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan]

[国内連絡先]

Phone : +81-3-5632-9617

Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Mailing address: Box 57, S-911 21 Vannas, Sweden
Phone/Fax: +46(0)935-33033/+46(0)935-33044

www.cosmobio.co.jp