

Product information

| | |
|------------------------|---|
| Product number: | AS05 061A |
| Lot Number: | 0603 |
| Clonality: | Polyclonal |
| Raised in: | Rabbit |
| Purity: | Affinity purified IgG on the antigen column |
| Quantity: | 200 µg |
| Antibody form: | Lyophilized In PBS pH 7.4 |

Please, add 200 µl of sterile dist. water for reconstitution of antibodies.

Storage instructions: 4°C temporary storage; -20°C or -80°C long term storage (years)

Please, avoid freezing and thawing of antibodies. Make aliquots instead.

Please, remember to spin tubes briefly prior to opening them to avoid any losses that might occur from liquid or lyophilized material adhering to the cap or sides of the tubes.

Background

Heat-shock protein 70 (Hsp70) is the major stress-inducible protein in vertebrates and highly conserved throughout evolution. It plays a role as a molecular chaperone and is important for allowing cells to cope with acute stressor insult, especially those affecting the protein machinery.

Immunogen: Synthetic peptide conjugated to KLH. Peptide sequence specific to C-terminal part of salmonid Hsp70. The target is not present in any other Hsp70 proteins.

Application information

Western blot: 1: 1000 needs to be optimised depending upon the system used.

Immunoprecipitation: 5 µg of antibodies have been used

Reactivity: The antibody is very specific for salmonid Hsp70 (Tested on salmon and brooktrout). It does not cross-react with Hsc70. It does not detect Hsp70 from other species.

Antibodies are intended for the research use only not for diagnostic or therapeutic use.

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN
URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan]

[国内連絡先]

Phone : +81-3-5632-9617

Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619