

Product information

Product number: AS05 079
Lot Number: 0601
Clonality: Polyclonal
Raised in: Hen
Purity: Total IgY
in PBS pH 7.4
Concentration: 28mg/ml
Antibody form: Liquid

Storage instructions: 4°C
*Preparation contains sodium azide
(0.02 %) as inhibitor of bacterial
growth.*

DO NOT FREEZE!

*Please, remember to spin tubes briefly prior to
opening them to avoid any losses that might occur
from liquid or lyophilized material adhering to the
cap or sides of the tubes.*

Related product:
Stanniocalcin-1 AS05 095
Suitable for a western blot technique.

Background


Stanniocalcin-1 is a homodimeric phosphoglycoprotein expressed in a wide variety of tissues, including kidney, prostate, thyroid, bone, and ovary. It is considered to have roles in many physiological processes, including bone development, reproduction, wound healing, angiogenesis, and modulation of inflammatory response.

Immunogen: Synthetic peptide conjugated to KLH.

Application information

Immunohistochemistry

5µg/ml for formalin-fixed and paraffin embedded tissue.
Antigen retrieval in 10 mM citric acid for 4x5 min at 620W in a microwave oven.



*Immunohistochemistry on normal human brain cortex.
The large neurons are strongly reactive.
a)abx20 b)non-immune controlx20
c)abx40 d)controlx40.*

Antibodies are intended for the research use only not for diagnostic or therapeutic use.

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

Mailing address: Box 57, S-911 21 Vannas, Sweden
Phone/Fax: +46(0)935-33033/+46(0)935-33044

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN
URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan]

[国内連絡先]

Phone : +81-3-5632-9617
FAX : +81-3-5632-9618

Phone : +81-3-5632-9610
FAX : +81-3-5632-9619