



COSMO BIO CO., LTD.
Inspiration for Life Science

Agrisera 
Antibodies for plant sciences

Anti-GnRHa AS05 088

Product information

Product number: AS05 088
Lot Number: 0602
Clonality: Polyclonal
Raised in: Hen
Purity: Total IgY
Quantity: 200 µl
Antibody form: Liquid
In PBS pH7.4

Storage instructions: 4°C
Preparation contains sodium azide (0.02 %) as inhibitor of bacterial growth.

DO NOT FREEZE!

Please, remember to spin tubes briefly prior to opening them to avoid any losses that might occur from liquid or lyophilized material adhering to the cap or sides of the tubes.

Background

Gonadotropin-releasing hormone (GnRH) is considered to stimulate LH secretion from the pituitary. This antibody is raised against an GnRH analog (GnRHa).

Immunogen: GnRHa (gonadotropin-releasing hormone agonist) conjugated to BSA

Application information

Potent GnRH agonist. Induce stimulation of LH secretion from the pituitary glands spawning activity in fish.

ELISA: 1 : 5 000

In this immunoassay the anti-GnRHa and sample/standard are added simultaneously within GnRHa coated wells. After incubation, the anti-GnRHa bound to the coated GnRHa is detected with HRP-labelled secondary antibodies and color development achieved by addition of TMB reagent.

Reactivity: N/A

Antibodies are intended for the research use only not for diagnostic or therapeutic use.

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN
URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan]

[国内連絡先]

Phone : +81-3-5632-9617

Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Mailing address: Box 57, S-911 21 Vannas, Sweden
Phone/Fax: +46(0)935-33033/+46(0)935-33044

www.cosmobio.co.jp