

**Product:** Anti-HSC70 (heat shock cognate protein 70)**Product no:** AS07 233

---

**Product Information****Antibody clonality:** Polyclonal**Raised in:** Rabbit**Purity:** Serum**Quantity:** 100 µl**Antibody form:** Lyophilized. For reconstitution please add 100 µl of sterile water.

Please, remember to spin tubes briefly prior to opening them to avoid any losses that might occur from liquid or lyophilized material adhering to the cap or sides of the tubes.

**Storage instructions:** -20°C or -80°C long term storage (years). Please, avoid freezing and thawing of antibodies. Make aliquots instead.

**Background**

The heat shock cognate protein 70 (HSC70) is a constitutively expressed house keeping protein which can act as a molecular chaperone towards other proteins.

**Immunogen:** HSC70 purified from wheat germ.**Application information:****Western Blot:** 1: 50 000 (with Alkaline Phosphatase system)**MW:** 70 kDa

*Antibodies are intended for the research use only not for diagnostic or therapeutic use.*

---

Product support: [inquiry@agrisera.com](mailto:inquiry@agrisera.com), <http://www.agrisera.com>