



Anti-Ppa2 (*S. pombe*) antibody, rabbit serum

BACKGROUND

Schizosaccharomyces pombe **Ppa2** is a type 2A-like serine/threonine-protein phosphatase catalytic subunit whose polypeptide sequence has ~80% identity to those of mammalian type 2A phosphatases. **Ppa2** determines the sensitivity to okadaic acid, which is an inhibitor of protein serine/threonine phosphatases. The loss of the **ppa2** gene causes cells to be hypersensitive to the okadaic acid. **Ppa2** plays important roles in cell cycle control. It may be involved in controlling the entry into mitosis, possibly acting as an inhibitor (ref.1). **Ppa2** is abundant in the cytoplasm, in contrast to the type 1-like phosphatase Dis2, which is enriched in the nucleus. Thus **Ppa2** may perform major functions outside the nucleus.

Product type	Primary antibodies
Host	Rabbit
Source	serum
Form	Liquid
	Rabbit antiserum added with 0.05 % sodium azide
Volume	100 µl
Concentration	
Specificity	Ppa2
Antigen	Recombinant C-terminal polypeptide (26kDa) of <i>S. pombe</i> Ppa2 (Ref. 1)
Isotype	

Application notes Immunoblotting, Immunofluorescence microscopy, Immunoprecipitation

Recommended use

Recommended dilutions

Immunoblotting (dilution: 1/1000)

Optimal dilutions/concentrations should be determined by the end user.

Data Link: Swiss-Prot [P23636](#)

Staining Pattern

Cross reactivity The antibody recognized both Ppa1 and Ppa2 polypeptides in *S. pombe* because of their high amino acid similarity (~80% identity) (Fig.1 and ref. 1).

Storage -20°C (long period, -70°C)

References 1) Kinoshita N *et al* "Negative regulation of mitosis by the fission yeast protein phosphatase ppa2." *Genes Dev* 7: 1059-1071 (1993) PMID: [8389306](#)

2) Kinoshita K *et al* "The regulatory subunits of fission yeast protein phosphatase 2A (PP2A) affect cell morphogenesis, cell wall synthesis and cytokinesis." *Genes Cell* 1:29-45 (1996) PMID: [9078365](#)




Fig.1 Identification of Ppa1 and Ppa2 proteins. An immunoblot with anti-ppa2 antibody is shown (ref.1).

lane 1: Wild-type *S. pombe*

lane 2: *ppa1*

lane 3: *ppa2*

lane 4: Wild-type carrying a multicopy plasmid with *ppa1* gene

lane 5: Wild-type carrying a multicopy plasmid with *ppa2* gene

lane 6: Wild-type carrying a multicopy plasmid with ADH promoter ligated with *ppa2* gene

The positions of *ppa1* (36 kDa) and *ppa2* (39 kDa) polypeptide bands are indicated.




Fig.2 Cellular location of Ppa1 and Ppa2 (ref.1).

Indirect immunofluorescence microscopy of *Δppa2* deletion, wild type (WT), and wild-type overexpressing *ppa2* (*nmt-ppa2*, 18hr) was done, using anti-ppa2 antibody (right);

The same cells stained by DAPI are also shown (left).

Immunofluorescence was hardly detected in *Δppa2* cells, whereas cytoplasmic immunofluorescence was abundant in wild-type cells. Wild-type cells carrying *nmt-ppa2* plasmid overexpress Ppa2 protein in the absence of thiamine for 18 hr. Immunofluorescence was enhanced further in the cytoplasm, often accumulated at the nuclear periphery or within restricted domains. The deformation of chromosomal DNA was also visible in overexpressed cells.

Bar, 10μm.

For research use only. Not for clinical diagnosis.

Manufactured by BioAcademia, Inc.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619