

## Anti-Vero Toxin1 (*E. coli*) / Shiga Toxin (*S. dysenteriae*) antibody, rabbit serum

64-025 100 µl

Vero toxin 1 (VT1) is produced by Vero toxin1-producing *E. coli* (VTEC) and has lethal activity to Vero cells. The primary structure of VT1 is identical or nearly identical to Shiga toxin (Stx) produced by *Shigella dysenteriae* serotype 1 and also called Slt 1 (Shiga-like toxin 1). VT1 is composed from one A subunit and five B subunits. Some *E. coli* strains produce both Slt1 and Slt2, and they share sequence identity of 55 %, but they are immunologically distinct.

To express the activity of VT/Stx, interaction with specific receptor Gb3 is indispensable. VT/Stx removes the 4324th adenine of 28S RNA of ribosome, inhibits protein synthesis and causes cell death. After invasion into cell subunit A is cut by furin to give A1 and A2. A1 is a catalytic fragment, and A2 is required for holo-enzyme formation by combining subunit B.

**Applications:** 1) Western blotting (2,000 fold dilution) (Fig. 1)

2) Immunoprecipitation (Fig. 2) 3) ELISA Other applications have not been tested.

**Immunogen:** Initial immunization by VT1 toxoid and booster by VT1 toxin.

**Reactivity:** VT1 of *E. coli* VTEC strain and Shiga toxin of *Shigella dysenteriae* 1.


**Form:** Rabbit antiserum added with 0.09% sodium azide.

**Storage:** Sent at 4°C. Upon arrival, spin-down, aliquot and store at -20°C

**Data link:** GenBank [M16625](#) Shiga-like toxin I subunit A and subunit B


UniProtKB/Swiss-Prot [Q9FBI2](#) Shiga toxin subunit A

UniProtKB/Swiss-Prot [Q7BQ98](#) Shiga toxin subunit B



**Fig.1.** Detection of VT1 by western blotting with anti-VT1 antibody.

1. SDS-PAGE of culture medium of VTEC, stained with CBB
2. SDS-PAGE of crude extracts of VTEC cells, stained with CBB
3. Western blotting of culture medium of VTEC
4. Western blotting of crude extracts of VTEC cells. Ani-VT1 antibody was used at 1/2,000 dilution



**Fig. 2.** Immunoprecipitation of VT1 from culture medium of VTEC with anti-VT1 antibody. Arrows shows subunit A and subunit B of VT1. Heavy chain and Light chain indicate those of IgG.

