

Anti-GRWD1 antibody, rabbit polyclonal, affinity-purified

70-130 100 ug

Key words: WD repeat, KIAA1942, WDR28, METTL18, Phosphorylation

Description: Glutamate-rich WD repeat-containing protein 1 (GRWD1) consists of 446 amino acids with molecular mass of 49.4 kDa. It has been found as a protein which interacts with METTL18 and CDT1 proteins. It has been implicated in regulation of DNA replication and/or in ribosome biogenesis.

Applications:

1. Western blotting (1/1,000~1/3,000 dilution)
2. Immunoprecipitation (Assay dependent)
3. Immunofluorescence staining / Immunochemistry (1/100~1/1,000 dilution)

Immunogen: Purified GST-GRWD1 (human, full-length) expressed in E. coli.

Reactivity: Human, mouse and rat. Other species have not been tested.

Purification: The antiserum was first adsorbed with GST-agarose column and then the pass-through fraction was affinity-purified with GST-GRWD1 agarose column.

Product: 1 mg/ml in PBS and 50% glycerol. Filter-sterilized. Carrier protein and azide free.

Storage: Sent at 4°C and upon arrival, spin-down and store at -20°C.

Database link: [uniprot/Q9BQ67](https://www.uniprot.org/entry/Q9BQ67) (GRWD1_HUMAN)

[uniprot/Q810D6](https://www.uniprot.org/entry/Q810D6) (GRWD1_MOUSE)

Reference: This product has been described and used in the following publication.

Sugimoto N. et al. (2008) Identification of novel human Cdt1-binding proteins by a proteomics approach: proteolytic regulation by APC/CCdh1. [Mol Biol Cell](https://doi.org/10.1007/s12013-008-9007-2). 19(3):1007-21.

WB, Free article.




Fig.1. Identification of GRWD1 proteins in whole cell lysates by western blotting with anti-GRWD1 antibody

Whole cell lysates of HeLa cells untreated (-) and treated (+) with DNA damaging agent, adriamycin (Ad), and NIH3T3 cells were analyzed by western blotting with anti-GRWD1 antibody at 1/1,000 dilution. The samples were 10 μ g. Second antibody was HRP-conjugated goat anti-rabbit IgG used at 1/5,000 dilution. The revelation of multiple bands indicates post-translational modification such as phosphorylation. The level of GRWD1 in the cell was not affected by DNA-damaging treatment. The identity of an additional band at ~85 kDa position other than the GRWD1 band in NIH-3T3 cell lysate is not known. The GRWD1 proteins were identified at a position higher (~55 kDa) than expected from the molecular mass of GRWD1 indicated from cDNA sequence (49.4 kDa).




Fig.2. Immunofluorescence staining of GRWD1 protein in HeLa cells.

HeLa cells were fixed in 4% paraformaldehyde overnight and permeabilized in 0.25% TritonX 100 in PBS for 10 min. Anti-GRWD1 antibody was used at 1/1,000 dilution. As second antibody, goat anti-rabbit IgG conjugated with Alex488 was used at 1/5,000 dilution. As a signal enhancer, Can Get Signal Immunostain B (Toyobo, Osaka) was used according to the protocol of the supplier. Nuclei were stained with DAPI. GRWD1 protein is localized in nuclei.