



## Anti-Dnmt1 (1-248)

### BACKGROUND

Dnmt1 (DNA (cytosine-5-)-methyltransferase 1) has a role in the establishment and regulation of tissue-specific patterns of methylated cytosine residues (epigenetics). Hypermethylation in promoter regions are associated with repression of the genes. Aberrant methylation patterns are associated with certain human tumors and developmental abnormalities. This antibody was raised by Prof. S. Tajima of Osaka Univ. and used in Ref.3

<b>Product type</b>	Primary antibodies
<b>Host</b>	Rabbit
<b>Source</b>	Serum
<b>Form</b>	Liquid 1 mg/ml in PBS, 50% glycerol, 0.05% sodium azide (and trace of ammonium sulfate) Affinity-purified with the recombinant Dnmt1
<b>Volume</b>	50ug
<b>Concentration</b>	1 mg/ml
<b>Specificity</b>	Reacts with Dnmt1 (1-248)
<b>Antigen</b>	Highly purified recombinant mouse Dnmt1 (amino acids 1-248)
<b>Isotype</b>	

**Application notes** WB, IP, IF Other applications were not tested

### Recommended use

Immunoprecipitation: Especially good for native form of Dnmt1 and useful for analysis of the protein complex containing Dnmt1

### Recommended dilutions

Western Blotting: 0.2~1 ug/ml

Optimal dilutions/concentrations should be determined by the end user.


### Staining Pattern

**Cross reactivity** Mouse and human. Not tested with other species

**Storage** -20 °C (long period, -80°C)

**References** This product was used in Ref. 3 and 4

- 1) Di Croce L, et al. (2002). Methyltransferase recruitment and DNA hypermethylation of target promoters by an oncogenic transcription factor. *Science* 295: 1079-82.
- 2) Rhee I, et al. (2002). DNMT1 and DNMT3b cooperate to silence genes in human cancer cells. *Nature* 416: 552-6.
- 3) Sharif K, Muto M, Takebayashi S, et al.(2007) The SRA protein Np95 mediates epigenetic inheritance by recruiting Dnmt1 to methylated DNA. *Nature* 450:908-12.



**Figure. Use of anti-Dnmt1 (1-248) antibody for western blotting and immuno-fluorescence staining**

- a, Subnuclear localization of Dnmt1 (green), Np95 (red) and PcnA (blue) in mouse embryonic stem cells (E14) during cell-cycle progression (bottom). Merged images for Np95 and Dnmt1 are shown. Profile of DNA content in exponentially growing ECSs is shown (top)
- b, Association of Np95 with Dnmt1 and PcnA in HeLa cell nuclear extracts. Human NP95 was tagged by biotin-binding domain and stably expressed in HeLa cells together with the *E. coli* BirA biotin ligase gene. Biotylated NP95 was captured by streptavidin beads. The captured proteins were eluted from the beads and analyzed by western blotting using antibodies against Dnmt1 and PcnA.
- c, Catalytic activity of NP95 associated Dnmt1 in HeLa cells. The presence of Dnmt1 in the NP95 complexes was confirmed by western blot and the DNA methyltransferase activity of the extracts was measured (The data were kindly provided by Dr. M. Muto of Riken Research Center for Allergy and Immunology)

**Related products:**

# 70-203 Anti-Dnmt1 (1037-1086) antibody, affinity-purified (rabbit polyclonal)

Especially good for immunoprecipitation of denatured Dnmt1 and useful for ChIP assays Cross-reacts with human, mouse and xenopus Dnmt1

*For research use only. Not for clinical diagnosis.*

Manufactured by BioAcademia Inc.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619