



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. BAM-70-415-EX

Anti-hnRNP-U / SAF-A antibody, rabbit serum

BACKGROUND

Heterogeneous nuclear ribonucleoprotein U (hnRNP-U), also known as **scaffold attachment factor A, SAF-A** is a nuclear matrix-associated protein that interacts with chromosomal DNA. **hnRNP-U** specifically binds to scaffold/matrix attachment region of DNA and could thus be involved in higher order chromatin structure. **hnRNP-U** is also a RNA binding protein and forms complexes with heterogeneous nuclear RNA (hnRNA) and plays an important role in pre-mRNA processing and transport.

HnRNP-U is reported to interact with necdin, a growth suppressor that is expressed in terminally differentiated neurons and skeletal muscle cells. It has been shown that **hnRNP-U** recruits necdin to the nuclear matrix where they form a stable complex. It is suggested that necdin suppresses cell proliferation through its interaction with **hnRNP-U** in the specific subnuclear structure (ref.2).

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
	Antiserum added with 0.05% sodium azide.
Volume	100 µl
Concentration	
Specificity	hnRNP-U / SAF-A
Antigen	Recombinant MBT-fused mouse hnRNP-U (aa 614-800)
Isotype	Mouse IgG2b kappa

Application notes WB, Immunocytochemistry, IP

Recommended use

Recommended dilutions

Western blotting (dilution: 1/3,000-1/1,000)
Immunocytochemistry (dilution: 1/1,000-1/500)

Optimal dilutions/concentrations should be determined by the end user.

Data Link: Swiss-Prot [Q8VEK3](#) (mouse), [Q00839](#) (human)

Cross reactivity

Reacts with mouse and rat, and predicted to react with human from the amino acid sequence homology.

Storage

-20°C

References

(This antibody is used in ref.2.)

- 1) Kiledjian M and Dreyfuss G (1992) "Primary structure and binding activity of the hnRNP U protein: binding RNA through RGG box." *EMBO J* 11: 2655-2664 PMID: [1628625](#)
- 2) Taniura H and Yoshikawa K (2002) "Necdin interacts with the ribonucleoprotein hnRNP U in the nuclear matrix." *J Cell Biochem* 84:545-555 PMID: [11813259](#)

Related Products

74-104 Anti-APP (C-terminal) antibody, rabbit serum




Fig.1 Immunoblotting of hnRNP-U with this antibody (ref.2).

Specificity of anti-hnRNP-U antibody, HUT.

Cell lysates were prepared from SAOS-2 cells transfected with pRc/CMV vectors (pRc) or pRc/CMV vectors expressing Myc-tagged hnRNP-U (Myc-UF). Exogenous Myc-tagged hnRNP-U (Myc-U) and endogenous hnRNP-U (U) proteins were detected by immunoblotting with anti-Myc antibody (αMyc) or HUT (αU).

This antibody recognized exogenous Myc-tagged hnRNP-U and endogenous ~120 kDa hnRNP-U proteins in SAOS-2 cells.




Fig.2 Immunocytochemistry using this antibody, HUT (ref.2)

Mouse P19 neurons were labeled with anti-necladin antibody (Ndn) (a) or with HUT for hnRNP-U (U) (c) in combination with anti-neuronal marker, MAP2, antibody for MAP2 (b, d). The nuclear matrix was prepared in situ and labeled for necladin (Ndn) (e), hnRNP-U (U) (f), and a nuclear matrix marker, lamin B (g).

Both necladin and hnRNP-U were localized to the nuclei of differentiated neurons, which express the neuronal marker MAP2 (a-d). Necladin was also distributed in the neuronal cytoplasm (a).

The immunocytochemical analysis of in situ extracted nuclear matrix revealed that both necladin and hnRNP-U were concentrated in intranuclear speckles throughout the nucleoplasm (e, f). Lamin B, a nuclear matrix marker, was localized to the nuclear lamina (g). These results suggest that both necladin and hnRNP-U are associated with the nuclear matrix of neurons.

For research use only. Not for clinical diagnosis.

Manufactured by BioAcademia, Inc.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp/index_e.asp e-mail: export@cosmobio.co.jp
[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610
FAX : +81-3-5632-9618 FAX : +81-3-5632-9619