

Anti-Nucleobindin 2 / NEFA / Nesfatin precursor antibody, rabbit polyclonal (NET1)

73-109 100 μ g

Storage: Ship at 4°C and store at -20°C.

Reactivity: Reacts with mouse rat and human nucleobindin 2. (Specificity validated in Ref)

Immunogen: Recombinant GST-fused mouse nucleobindin 2 (aa 26-420)

Applications:

1. Western blotting (1/1,000-1/3,000).
2. Immuno-precipitation. (1/200)
3. Immuno-cytochemistry (1/300-1/1,000)
4. Immuno-histochemistry (1/300-1/1,000)
5. Immuno-electron microscopy (assay dependent)
6. Immuno-affinity chromatography (assay dependent)

Form: Protein A-affinity .purified IgG. 1 mg/ml in PBS, 50% glycerol. Filter-sterilized. No additive.

Background: **Nucleobindin 2 (NUCB2)**, also known as **NEFA** or **Nesfatin precursor**, is a ubiquitously expressed EF-hand Ca^{2+} binding protein that is implicated in various physiological processes.


Nucleobindin 2 interacts with the postmitotic growth suppressor neclin in neurons. Both neclin and **nucleobindin 2** are expressed in differentiated neurons and skeletal muscles and these proteins are likely to be involved in the regulation of survival and death of postmitotic cells by controlling Ca^{2+} homeostasis

Data Link: Swiss-Prot [P81117](#) (mouse), [Q9JI85](#) (rat), [P80303](#) (human)

Reference: This antibody was described and used in the following publication.

Taniguchi N *et al* (2000) "The postmitotic growth suppressor neclin interacts with a calcium-binding protein (NEFA) in neuronal cytoplasm." *J Biol Chem* **275**: 31674-31681 PMID: [10915798](#) **WB, IP,, IF,**

IHC, Immuno-affinity chromatography






Fig.2 Expression of Nucleobindin-2 in various tissues as examined by western blotting.

Distribution of nucleobindin 2 in neonatal mouse organs. Homogenates of various organs from P0 mouse were separated by 10% SDS-PAGE and immunoblotted with this antibody. The antibody was used at 1/1,000 dilution.




Fig.3. Immunocytochemistry for endogenous nucleobindin 2

Cells were stained with this antibody by the avidin-biotin-peroxidase complex method.

Left panel; undifferentiated murine embryonal carcinoma P19 cells (UD). Right panels, enriched post-mitotic neurons (PN). Nucleobindin 2 was localized to the cytoplasm near the nucleus in undifferentiated P19 cells, and its immunoreactivity in the cytoplasm was increased when P19 cells were induced to differentiate into neurons.




Fig.4 Immunohistochemistry for nucleobindin 2 in neonatal mouse brain with this antibody

Frozen brain sections from neonatal mouse were stained with this antibody by the avidin-biotin peroxidase complex method. A-C, cerebral cortex (parietal lobe). At higher magnification (B, C), fine granular immunoreactive materials are observed at both neuronal dendrites (arrows) and perikarya (asterisks) in the layer IV (arrowhead B in A) and subplate (arrowhead C in A) of the cerebral cortex. Scale bars, 100 um (A) and 10um (B and C).