

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-AK2-286

Anti AK2

Product type Primary antibodies

Immunogen AK2 (NP_001616.1, 1 a.a. ~ 239 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MAPSVPAAEPEYPKGIRAVLLGPPGAGKGTQAPRLAENFCVCHLATGDMLRAMVASGSEL
GKKLKATMDAGKLVSDVMVELIEKNLETPLCKNGFLLDGFPRTVRQAEMLDDLMEKRKE
KLDSVIEFSIPDSLLIRITGRLIHPKSGRSYHEEFNPPKEPMKDDITGEPLIRRSDDNE
KALKIRLQAYHTQTTPLIEYYRKRGIHSAIDASQTPDVVFASILAAFSKATCKDLVMFI

Gene Alias ADK2

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~26.50 kDa using AK2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of AK2 expression in transfected 293T cell line by AK2 rabbit polyclonal antibody.

Lane 1: AK2 transfected lysate(26.50 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of AK2 expression in human kidney (Predicted Mw: 26.3470 kDa).

Western Blotting (Tissue lysate)



AK2 rabbit polyclonal antibody. Western Blot analysis of AK2 expression in human kidney.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619