



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-ALDH7A1-225

Anti ALDH7A1

Product type Primary antibodies

Immunogen ALDH7A1 (NP_001173.1, 1 a.a. ~ 511 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MSTLLINQPQYAWLKELGLREENEGVYNGSWGGRGEVITYCPANNEPIARVRQASVADY
EETVKKAREAWKIWADIPAPKRGEIVRQIGDALREKIQVLGSLVSLEMGKILVEGVGEVQ
EYVDICDYAVGLSRMIGGPILPSERSGHALIEQWNPVGLVGIITAFNFPVAVYGWNNIAIA
MICGNVCLWKGAPTTSLISVAVTKIIAKVLEDNKLPGAICSLTCGGADIGTAMAKDERVN
LLSFTGSTQVGKQVGLMVQERFGRSLLELGGNNIAIFEDADLSLVVPSALFAAVGTAGQ
RCTTARRLFIHESIHDDEVVNRLKKAYAQIRVGNPWPDPNVLYGPLHTKQAVSMFLGAVEEA
KKEGGTVVYGGKVMMDRPGNYVEPTIVTGLGHDASIAHTETFAPILYVFKFKNEEEVFAWN
NEVKQGLSSSIFTKDLGRIFRWLGPKGSDCGIVNVNIPTSGAEIGGAFGGEKHTGGGRES
GSDAWKQYMRRSTCTINYSKDLPLAQGIKFQ

Gene Alias ATQ1, EPD, FLJ11738, FLJ92814, PDE

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~55.40 kDa using ALDH7A1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of ALDH7A1 expression in transfected 293T cell line by ALDH7A1 rabbit polyclonal antibody.

Lane 1: ALDH7A1 transfected lysate(55.40 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619