



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-ALG2-245

Anti ALG2


Product type	Primary antibodies
Immunogen	ALG2 (NP_149078.1, 1 a.a. ~ 416 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence
MAEEQGRERDSVPKPSVLFLHPDLGVGGAERLVLDAALALQARGCSVKIWTAHYDPGHCF
AESRELPVRCAGDWLPRGLGWGGRGAAVCAYVRMVFLALYVFLADEEFDVVVCDQVSAC
IPVFRLARRRKILFYCHFPDLLLTKRDSFLKRLYRAPIDWIEEYTTGMADCILVNSQFT
AAVFKETFKSLSHIDPDVLYPSLNVTSFDSVVPEKLDLVPKGKKFLLLSINRYERKKNL
TLALEALVQLRGRLLTSQDWERVHLIVAGGYDERVLENVEHYQELKMMVQQSDLGQYVTFLL
RSFSDKQKISLLHSCVLYTPSNEHFGIVPLEAMYMQCPVIAVNSGGPLESIDHSVTGF
LCEPDPVHFSEAIEKFIREFSLKATMGLAGRARVKEKFSPEAFTEQLYRYVTKLLV

Gene Alias CDGLI, FLJ14511, hALPG2

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~47.10 kDa using ALG2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of ALG2 expression in transfected 293T cell line by ALG2 rabbit polyclonal antibody.

Lane 1: ALG2 transfected lysate(47.10 KDa).

Lane 2: Non-transfected lysate.

.Antibody dilution: 1:1000. Western Blot analysis of ALG2 expression in A-431 (Predicted Mw: 47.0920 kDa).


Western Blotting (Cell lysate)



ALG2 rabbit polyclonal antibody. Western Blot analysis of ALG2 expression in A-431.

Antibody dilution: 1:1000. Western Blot analysis of ALG2 expression in human kidney (Predicted Mw: 47.0920 kDa).

Western Blotting (Cell lysate)



ALG2 rabbit polyclonal antibody. Western Blot analysis of ALG2 expression in A-431.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619