

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-AMY1C-467

Anti AMY1C

Product type Primary antibodies

Immunogen AMY1C (AAI56581.1, 1 a.a. ~ 511 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MKLFWLLFTIGFCWAQYSSNTQQGRTSIVHLFEWRWVDIALECERYLAPKGFGGVQVSP
NENVAIHNPFRPWWERYQPVS YKLCTRSGNEDEFNMVTRCANNVGVRIYVDAVINHMCGN
AVSAGTSSTCGSYFNPGSRDFPAVPYSGWDFNDGKCKTGSGDIENYNDATQVRDCRLSGL
LDLALGKDYVRSKIAEYMNHLIDIGVAGFRIDASKHMPGDIKAILDKLHNLNSNWFPEG
SKPFIYQEVIDLGGEPIKSSDYFGNGRVTEFKYGAKLGTVIRKWNGEKMSYLKNWGEGWG
FMPSDRALVFVDNHDNQRGHGAGGASILTFWDARLYKMAVGFM LAHPYGFTRVMSSYRWP
RYFENGKDVNDWVGPPNDNGVTKEVTINPDTCGNDWVCEHRWRQIRNMVNFRNVVDGQP
FTNWDNGSNQVAFGRGNRGFIVFNDDWTFSLTLQTGLPAGTYCDVISGDKINGNCTGI
KIYVSDDGKAHFSISNSAEDPFIAIHAESKL

Gene Alias AMY1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~56.21 kDa using AMY1C transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of AMY1C expression in transfected 293T cell line by AMY1C rabbit polyclonal antibody.

Lane 1: AMY1C transfected lysate(56.21 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of AMY1C expression in HeLa (Predicted Mw: 57.7680 kDa).


Western Blotting (Cell lysate)



AMY1C rabbit polyclonal antibody. Western Blot analysis of AMY1C expression in HeLa.

Antibody dilution: 1:1000. Western Blot analysis of AMY1C expression in human pancreas (Predicted Mw: 57.7680 kDa).

Western Blotting (Tissue lysate)



AMY1C rabbit polyclonal antibody. Western Blot analysis of AMY1C expression in human pancreas.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619