

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-AOC2-265


Anti AOC2

Product type	Primary antibodies
Immunogen	AOC2 (NP_033720.2, 1 a.a. ~ 756 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MHLKIVLAFLALSLITIFALAYVLLTSPGGSSQPPHCPSVSHRAQPWPHPGQSQLFADLS REELTAVMRFLTQRLGPGLVDAQAQPSDNCIFSVELQLPPKAAALAHLDRGSPPPAREA LAIVLFGGQPQPNVSELVVGPLPHPSYMRDVTVERHGGPLPYHRRPVLRAEFTQMWRHLK EVELPKAIFLSSTFNYNGSTLAHVHATPRGLRSGDRATWMALYHNISGVGLFLHPVGLE LLLDHRALDPAHWTVQQVFYLGHYADLGQLEREFKSGRLEVVRVPLPPPNGASSLRSRN SPGPLPPLQFSPQGSQYSVQGNLVVSSLWSFTFGHGVFSGLRIFDVRVFQGERIAYEVSQ ECVSIYGADSPKTMTRYLDSSFGLGRNSRGLVRGVDCPYQATMVDIHILVGKGAVQLLP GAVCVFEEAQGLPLRRHHNYLQNHFYGGGLASSALVRSVSSVGNVDYIWDVFLYPNGALE GRVHATGYINTAFLKGGEELLFGNRVGERVLGTVHTHAFHFKLDLDVAGLKNWVVAEDV VFKPVAAPWNPEHWLQRPQLTRQVLGKEDLTAFLGSPLPRYLALASNQTNAWGHQRGYR IQIHSPGLGIHPILESDMERALSWGRYQLVVTQRKEEESQSSSIYHQNDIWTPTVTFADFI NNETLLGEDLVAVVTASFLHIPHAEDIPNTVTLGNRVGFLLRPYNFFDEDPSIFSPGSVY FEKGQDAGLCSINPVACLPDLAACVPDLPPFSYHGF

Gene Alias DAO2, RAO

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~83.70 kDa using AOC2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of AOC2 expression in transfected 293T cell line by AOC2 rabbit polyclonal antibody.

Lane 1: AOC2 transfected lysate(83.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619