

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-APOBEC3C-090


Anti APOBEC3C

Product type	Primary antibodies
Immunogen	APOBEC3C (NP_055323.2, 1 a.a. ~ 190 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MNPQIRNPMKAMYPGTFYFQFKNLWEANDRNETWLCFTVEGIKRRSVVSWKTGVFRNQVD SETHCHAERCFLSWFCDDILSPNTKYQVTWYTSWSPCPDCAGEVAEFLARHSNVNLTIFT ARLYYFQYPCYQEGLRSLSQEGVAVEIMDYEDFKYCWENFVYNDNEPFKPKWGLKTNFRL LKRRRLRESLQ

Gene Alias APOBEC1L, ARDC2, ARDC4, ARP5, MGC19485, PBI, bK150C2.3

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~22.80 kDa using APOBEC3C transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of APOBEC3C expression in transfected 293T cell line by APOBEC3C rabbit polyclonal antibody.


Lane 1: APOBEC3C transfected lysate(22.80 kDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of APOBEC3C transfected lysate.

Immunoprecipitation



Immunoprecipitation of APOBEC3C transfected lysate using APOBEC3C rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with APOBEC3C purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619