

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-ARHGEF10-332


Anti ARHGEF10

Product type	Primary antibodies
Immunogen	ARHGEF10 (AAH40474.1, 1 a.a. ~ 380 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MENPEEAIYDDVPRENSDSEPDEMIYDDVENGDEGGNSSLEYGWSSEFESYEEQSDSEC KNGIPRSFLRSNHKKQMQLVKAAKDGTGDGLERTRAAVKRGRSFIRTKSLIAQDHRSSL EEEQNLFDVDCKHPEAILTPMPEGLSQQQVRRYILGSVVDSEKNYVDALKRILEQYK PLSEMEPKVLSEKRLKTVFYRVKEILQCHSLFQIALASRVSEWDSVEMIGVVFVASFSKS MVLDAYSEYVNNFSTAVAVLKKTCAKPAFLKQEQEASPDRTTLYSLMMKPIQRFPQ FILLQDMLKNTSKGHPDRLPLQMALTELETLAEKLNKRDADQRCEVKQIAKAINERY LNKVERGFLQLYSKIIFALC

Gene Alias DKFZp686H0726, GEF10, MGC131664

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~43.70 kDa using ARHGEF10 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of ARHGEF10 expression in transfected 293T cell line by ARHGEF10 rabbit polyclonal antibody.

Lane 1: ARHGEF10 transfected lysate(43.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619