

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-ARID3A-577

Anti ARID3A

Product type	Primary antibodies
Immunogen	ARID3A (AAH60828.1, 1 a.a. ~ 593 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MKLQAVMETLLQRQQRARQLEARQQLPPDPPAAPHGRARAAPDEDREPESARMQRAQMA ALAAMRAAAAGLGHPASPGGSEDGPPGSEEDAAAREGTPGSPGRGREGPGEEHFEDMASD EDMKPKWEEEEEMEEDLGEDEEEEEEDYEDEEEEEDEEGLGPPGPASLGTTALFPRKAQPP QAFRGDGVPRVLGGQERPGPGPAHPGGAAHVAPQLQPPDHGDWYEEQFKQLYELDGDPK RKEFLDDLFSFMQKRGTTPVNRIPIMAKQVLDLFMLYVLVTEKGGLVEVINKKLWREITKG LNLPTSITSAFTLRTQYMEYLYPYECEKRGLSNPNELQAIDSNRREGRRQSFGGSLFA YSPGGAHGMLSSPKLPVSSLGLAASTNGSSITPAPKIKKEEDSAIPITVPGRLPVSLAGH PVVAAQAAAVQAAAAQAQAAQAAALEQLREKLESAEPPEKKMALVADEQQRLMQRALQQ NFLAMAAQLPMSIRINSQASESRQDSAVNLTGTNGSNSISMSVEINGIMYTGVLFAQPPA PTPTSAPNKGSGGGGGSSSSNAGGRGGNTGTSGGQAGPAGLSTPSTSTSNNSLP
Gene Alias	BRIGHT, DRIL1, DRIL3, E2FBP1
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~63.00 kDa using ARID3A transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of ARID3A expression in transfected 293T cell line by ARID3A rabbit polyclonal antibody.

Lane 1: ARID3A transfected lysate(63.00 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619