

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-BBS4-410

Anti BBS4

Product Primary antibodies

type

Immunogen BBS4 (NP_149017.2, 1 a.a. ~ 519 a.a) full-length human protein..

Form Liquid

Protein A purified

Volume 100ul

Specificity Human

Host Rabbit


Immunogen Sequence MAEERVATRTQFPVSTESQKPRQKKAPEFPILEKQNWLIHLHYIRKDYEACKAVIKEQLQ
ETQGLCEYAIYVQALIFRLEGNIQESLELFQTCAVLSPQSADNLKQVARSLFLLGKHKAA
IEVYNEAAKLNQKDWEISHNLGVCYIYLKQFNKAQDQLHNALNLNRHDLTYIMLGKIHLL
EGDLDKAIEVYKKAVEFSPENTELLTTLGLLYLQLGIYQKAFEHLGNALTYDPTNYKAIL
AAGSMMQTHGDFDVALTKYRVVACAVPESPPLWNNIGMCFFGKKKYVAAISCLKRANYLA
PFDWKILYNLGLVHVTMQQYASAFHFLSAINFQPKMGELYMLLAVALTNLEDIENAKRA
YAEAVHLDKCNPLVNLNYAVLLYNQGEKKNALAQQEMEKVSLKDNSSLEFDSEMVEM
AQKLGAALQVGEALVWTKPVKDPKSKHQTSTSKPASFQQPLGSNQALGQAMSSAAAYRT
LPSGAGGTSQFTKPPSLPLEPEPAVESSPTETSEQIREK

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~58.30 kDa using BBS4 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of BBS4 expression in transfected 293T cell line by BBS4 rabbit polyclonal antibody.

Lane 1: BBS4 transfected lysate(58.30 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of BBS4 expression in mouse liver (Predicted Mw: 58.2820 kDa).

Western Blotting (Tissue lysate)




BBS4 rabbit polyclonal antibody. Western Blot analysis of BBS4 expression in mouse liver.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of BBS4 transfected lysate.

Immunoprecipitation



Immunoprecipitation of BBS4 transfected lysate using BBS4 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with BBS4 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619