



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CBLC-367

Anti CBLC

Product type	Primary antibodies
Immunogen	CBLC (AAH28915.1, 1 a.a. ~ 474 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MALAVAPWGRQWEEARALGRAVRMLQRLEEQCVDPRLSVSPPSLRDLLPRTAQLLREVAH SRRAAGGGGPGGGPGGSGDFLLIYLANLEAKSRQVAALLPPRGRRSANDELFRAGSRLRRQ LAKLAIIFSHMHAELHALFPGGKYCGHMYQLTKAPAHTFWRESCGARCVPWAEFESLLG TCHPVEPGCTALALRTTIDLTCSGHVSIFEFDVFTRLFQWPPTLLKNWQLLAVNHPGYMA FLTYDEVQERLQACRDKPGSYIFRPSCTRLGQWAIGYVSSDGSILQTIPANKPLSQVLLE GQKDGIFYLYPDGKTHNPDLTELGQAEPQQRIHVSEEQLQLYWAMDSTFELCKICAESNKD VKIEPCGHLLCSCCLAAWQHSDSQTCPFCRCCEIKGWEAVSIYQFHGQATAEDPGNSSDQE GRELELGQVPLSAPPLPPRPDLPPRKPRNAQPKVRLLKGNSPPAALGPQDPAPA


CBL-3, CBL-SL, RNF57

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~52.50 kDa using CBLC transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CBLC expression in transfected 293T cell line by CBLC rabbit polyclonal antibody.

Lane 1: CBLC transfected lysate(52.50 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619