

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*


**Catalog No. CNP-CCL18-040**

# Anti CCL18

|                           |                                                                                               |
|---------------------------|-----------------------------------------------------------------------------------------------|
| <b>Product type</b>       | Primary antibodies                                                                            |
| <b>Immunogen</b>          | CCL18 (NP_002979.1, 1 a.a. ~ 89 a.a) full-length human protein.                               |
| <b>Form</b>               | Liquid                                                                                        |
|                           | Protein A purified                                                                            |
| <b>Volume</b>             | 100ul                                                                                         |
| <b>Specificity</b>        | Human                                                                                         |
| <b>Host</b>               | Rabbit                                                                                        |
| <b>Immunogen Sequence</b> | MKGLAAALLVLVCTMALCSCAQVGTNKELCCLVYTSWQIPQKFIVDYSETSPQCPKPGVI<br>LLTKRGRQICADPNKKWVQKYISDLKLNA |
| <b>Gene Alias</b>         | AMAC-1, AMAC1, CKb7, DC-CK1, DCCK1, MIP-4, PARC, SCYA18                                       |

**Application notes** **WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~9.80 kDa using CCL18 transfected 293T lysate.

**Western Blotting (Transfected lysate)**



Western Blot analysis of CCL18 expression in transfected 293T cell line by CCL18 rabbit polyclonal antibody.


Lane 1: CCL18 transfected lysate(9.80 kDa).

Lane 2: Non-transfected lysate.

## IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of CCL18 transfected lysate.

### Immunoprecipitation



Immunoprecipitation of CCL18 transfected lysate using CCL18 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with CCL18 rabbit polyclonal antibody.

## Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619