

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-CCM2-170

Anti CCM2

Product type	Primary antibodies
Immunogen	CCM2 (NP_113631.1, 1 a.a. ~ 444 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MEEEGKKGKKPGIVSPFKRVFLKGEKSRDKKAHEKVTERRP LHTVVLSPERVEPDRLLS DYIEKEVKYLGQLTSIPGYLNPSSRTEILHFIDNAKRAHQ LPGHLTQEHDVLSLSAYNV KLAWRDGEDILRVPIHDIAAVSYVRDDAAHLVVLKTAQDPGISPSQSLCAESSRGLSAG SLSESAVGPVEACCLVILAAESKVAEEELCCLLGQVFQVVYTESTIDFLDRAIFDGASTP THHLSLHSDDSSTKVDIKETYVEASTFCFPESVDVGGASPHSKTISESELSASATELLQ DYMLTLRTKLSSQEIQQFAALLHEYRNGASIHEFCINLRQLYGDSRKFLLLGLRPFPIEK DSQHFENFLETIGVKDGRGIITDSFGRHRRALSTTSSSTTNGNRATGSSDDRSAPSEGDE WDRMISDISSDIEALGCSMDQDSA
Gene Alias	C7orf22, MGC4067, MGC4607, MGC74868, PP10187
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~48.80 kDa using CCM2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CCM2 expression in transfected 293T cell line by CCM2 rabbit polyclonal antibody.

Lane 1: CCM2 transfected lysate(48.80 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619