

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CHIA-088

Anti CHIA


Product type	Primary antibodies
Immunogen	CHIA (NP_068569.2, 1 a.a. ~ 368 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MVSTPENRQTFITSVIKFLRQYEFDGLDFDWEYPGSRGSPQDKHLFTVLVQEMREAFEQ EAKQINKPRLMVTAAVAAGISNIQSGYEIPQLSQYLDYIHVMYDLHGSWEGYTGENSPL YKYPTDTGSNAYLNVDYVMNYWKDNGAPAEKLIVGFPTYGHNFILSNPSNTGIGAPTSGA GPAGPYAKESGIWAYYEICTFLKNGATQGWDAPQEVYAYQGNVWVGYDNIKSFDIKAQW LKHNFEGGAMVWAILDDFTGTFCNQKGFPLISTLKKALGLQSASCTAPAQPIEPITAAP SGSGNGSGSSSGSGSGSGFCAVRANGLYPVANNRNAFWHCVNGVTYQQNCQAGLVFDT SCDCCNWA
---------------------------	--

Gene Alias	2200003E03Rik, AMCase, CHIT2, DKFZp313J1722, RP5-1125M8.1, TSA1902
-------------------	--

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~40.10 kDa using CHIA transfected 293T lysate.
--------------------------	--

Western Blotting (Transfected lysate)




Western Blot analysis of CHIA expression in transfected 293T cell line by CHIA rabbit polyclonal antibody.

Lane 1: CHIA transfected lysate(40.10 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of CHIA expression in mouse intestine (Predicted Mw: 52.2710 kDa).

Western Blotting (Tissue lysate)



CHIA rabbit polyclonal antibody. Western Blot analysis of CHIA expression in mouse intestine.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619