

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CKMT2-083

Anti CKMT2

Product type	Primary antibodies
Immunogen	CKMT2 (AAH29140.1, 1 a.a. ~ 419 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit


Immunogen Sequence

MASIFSKLLTGRNASLLFATMGTSVLTTGYLLNRQKVCAEVREQPRLFPPSADYPDLRKH
 NNCMAECLTPAIYAKLRNKVTPNGYTLDQCIQTGVDNPGHPFIKTVGMVAGDEESYEVFA
 DLFDPVIKLRHNGYDPRVMKHTTDL DASKITQQGFDEHYVLSSRVRTGRSIRGLSLPPAC
 TRAERREVENVAITALEGLKGDLAGRYYYKLSEMTEQDQQRLIDDHFLDKPVSPLLT CAG
 MARDWPDARGIWHNYDKTFLIWINEEDHTRVISMEEKGGNMKRVFERFCRGLKEVERLIQE
 RGWEFMWNERLG YILTCPSNLGTGLRAGVHVRIPKLSKDPRFSKILENLR LQKRGTTGGVD
 TAAVADVYDISNIDRIGRSEVELVQIVIDGVNYLVDCEKKLERGQDIKVPPLPQFGKK

Gene Alias SMTCK

Application notes **WB**
 Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~47.50 kDa using CKMT2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CKMT2 expression in transfected 293T cell line by CKMT2 rabbit polyclonal antibody.

Lane 1: CKMT2 transfected lysate(47.50 KDa).

Lane 2: Non-transfected lysate.

1

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619