

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CMPK-281

Anti CMPK

Product type Primary antibodies

Immunogen CMPK1 (AAH14961.1, 1 a.a. ~ 228 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MLSRCRSRLHVLGLSFLQTRRPILLCSPRLMKPLVVFVLGGPGAGKGTQCARIVEKYG
YTHLSAGELLRDERKNPDSQYGELIEKYIKEGKIVPEITISLLKREMDQTMAANAQKNK
FLIDGFPRNQDNLQGWNKTMGKADVSVLFFDCNNEICIERCLERGKSSGRSDDNRESL
EKRIQTYLQSTKPIIDLYEEMGKVKKIDASKSVDEVFDEVVQIFDKEG

Gene Alias CMK, CMPK, UMK, UMP-CMPK, UMPK

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~26.00 kDa using CMPK1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CMPK1 expression in transfected 293T cell line by CMPK1 rabbit polyclonal antibody.

Lane 1: CMPK1 transfected lysate(26.00 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619