

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CRLF3-369

Anti CRLF3

Product type	Primary antibodies
Immunogen	CRLF3 (AAH23567.2, 1 a.a. ~ 442 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MRGAMELEPELLLQEARENVEAAQSYRRELGHRLEGLREARRQIKESASQTRDVLKQHFN DLKGTGLGKLLDERLVTLLQEVDITIEQETIKPLDDCQKLIHGVNTAEDLVREGEIAMLGG VGEENEKLWSFTKKASHIQLDSLPEVPLLVDVPCLSAQLDDSIILNIVKDHIFKHGTVASR PPVQIEELIEKPGGIIVRWCKVDDDFTAQDYRLQFRKCTSNHFEDVYVGSETEFIVLHID PNVDYQFRVCARGDGRQEWSPWSVPQIGHSTLVPHEWTAGFEGYSLSSRRNIALRNDSES SGVLYSRAPTYFCGQTLTFRVETVGQPDRRDSIGVCAEKQDGYDSLQRDQAVCISTNGAV FVNGKEMTNQLPAVTSGSTVTFDIEAVTLGTTSNNEGGHFKLRVTISSNNREVVFDWLLD QSCGSLYFGCSFFYPGWKVLVF


CREME9, CYTOR4, FRWS, MGC20661

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~49.80 kDa using CRLF3 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of CRLF3 expression in transfected 293T cell line by CRLF3 rabbit polyclonal antibody.

Lane 1: CRLF3 transfected lysate(49.80 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of CRLF3 expression in mouse liver (Predicted Mw: 49.7660 kDa).

Western Blotting (Tissue lysate)




CRLF3 rabbit polyclonal antibody. Western Blot analysis of CRLF3 expression in mouse liver.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of CRLF3 transfected lysate.

Immunoprecipitation



Immunoprecipitation of CRLF3 transfected lysate using CRLF3 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with CRLF3 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619