

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CTBP1-478

Anti CTBP1

Product type Primary antibodies

Immunogen CTBP1 (NP_001012632.1, 1 a.a. ~ 429 a.a) full-length human protein

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MSGVRPPIMNGPLHPRPLVALLDGRDCTVEMPILKDVATVAFCDQAQSTQEIHEKVLNEAV
GALMYHTITLTREDLEKFKALRIIVRIGSGFDNIDIKSAGDLGIAVCNVPAASVEETADS
TLCHILNLYRRATWLHQALREGTRVQSVEQIREVASGAARIRGETLGIIGLGRVGQAAVAL
RAKAFGFNVLFYDPYLSGVERALGLQRVSTLQDLLFHSDCVTLHCGLNEHNNHLLINDFT
VKQMRQGAFLVNTARGGLVDEKALAQAQALKEGRIRGAALDVHESEPFSFSQGPKDAPNLI
CTPHAAWYSEQASIEMREEAAREIRRAITGRIPDSLKNCVNKDHLTAATHWASMDPAVVH
PELNGAAYRYPPGVGVAPTGIPAAVEGIVPSAMSLSHGLPPVAHPPHAPSPGQTVKPEA
DRDHASDQL

Gene Alias BARS, MGC104684

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~46.40 kDa using CTBP1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CTBP1 expression in transfected 293T cell line by CTBP1 rabbit polyclonal antibody.

Lane 1: CTBP1 transfected lysate(46.40 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619