

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CTBP2-573

Anti CTBP2

Product type Primary antibodies

Immunogen CTBP2 (AAH37900.1, 1 a.a. ~ 513 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MRPGLPGPTGLCAQTSSRGQKSVLKQKESCGIWQLYHFLSRKQEPRWEPCVSGSSSGDGA
VADLADELRGYPALCCTLPVHSYRSWAGIRPQIMNGPLHPRPLVALLDGRDCTVEMPILK
DLATVAFCDAAQSTQEIHEKVLNEAVGAMMYHTITLTREDLEKFKALRVIVRIGSGYDNVD
IKAAGELGIAVCNIPSAAVEETADSTICHILNLYRRNTWLYQALREGTRVQSVEQIREVA
SGAARIRGETLGLIGFGRTGQAVAVRAKAFGFSVIFYDPYLQDGIERSLGVQRVYTLQDL
LYQSDCVSLHCNLNEHNNHLLINDFTIKQMRQGAFLVNAARGGLVDEKALAQALKEGRIRR
AALDVHESEPFSAQGPLKDAPNLICTPHTAWYSEQASLEMREAAATEIRRAITGRIPES
LRNCVNKEFFVTSAPWSVIDQQAHPHPELNGATYRYPPGIVGVAPGGLPAAMEGIIPGGIP
VTHNLPTVAHPSQAPSPNQPTKHGDNREHPNEQ


Gene Alias

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~56.20 kDa using CTBP2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CTBP2 expression in transfected 293T cell line by CTBP2 rabbit polyclonal antibody.

Lane 1: CTBP2 transfected lysate(56.20 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619