

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-CTNND1-140

Anti CTNND1

Product type	Primary antibodies
Immunogen	CTNND1 (AAH10501.1, 1 a.a. ~ 610 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MIGEEVPSDQYYWAPLAQHERGSLASLDSLRKGGPPPPNWRQPELPEVIAMLGFRLLDAVK SNAAAYLQHLCYRNDKVKTDVRKLKGIPVLVGLLDHPKKEVHLGACGALKNISFGRDQDN KIAIKNCDGVPALVRLLRKARDMDLTEVITGTLWNLSHDSIKMEIVDHALHALTDEVII PHSGWEREPNEDCKPRHIEWESVLTNTAGCLRNVSSESEARRKLRECDGLVDALIFIVQ AEIGQKDSDSLKLVENCVCLLRNLSYQVHREIPQAERYQEAAPNVANNTGPHAASCFGAKK GKGKKPIEDPANDTVDFPKRTSPARGYELLFQPEVVRIYISLLKESKTPAILEASAGAIQ NLCAGRWTYGRYIRSALRQEKAISAIADLLTNEHERVVKAASGALRNLAVDARNKELIGK HAIPNLVKNLPGGQQNSSWNFSEDTVISILNTINEVIAENLEAAKKLRETQGIEKLVLIN KSGNRSEKEVRAAALVLQTIWGYKELRKPLEKEGWKKSDFFQVNLNNASRSQSSHSYDDST LPLIDRNQKSDKKPDREEIQMSNMGSNTKSLDNNYSTPNERGDHNKTLDRSGDLGDMPEPL KGTTPLMQKI
Gene Alias	CAS, CTNND, KIAA0384, P120CAS, P120CTN, p120
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~67.10 kDa using CTNND1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CTNND1 expression in transfected 293T cell line by CTNND1 rabbit polyclonal antibody.

Lane 1: CTNND1 transfected lysate(67.10 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619