

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CTPS2-400

Anti CTPS2

Product type	Primary antibodies
Immunogen	CTPS2 (NP_062831.3, 1 a.a. ~ 586 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MKYILVTGGVISGIGKGIASSIGTILKSCGLRVTAIKIDPYINIDAGTFSPYEHGEV FV LNDGGEVDLDLGN YERFLDINLYKDN NITTKIYQHVINKERRGDYLGKTVQVPHITDA VQEWVMNQAKVPVDGNKEEPQICVIELGGTIGDIEGMPFVEAFRQFQFKAKRENFCNIHV SLVPQLSATGEQKTKPTQNSVRALRGLGLSPDLIVCRSSTPIEMAVKEKISMFCNVNPEQ VICIHDVSSTYRVPVLL EEQSIVKYFKERLHLPIGDSASNLLFKWRNMADRYERLQKICS IALVGKYTKLRDCYASVFKALEHSALAINHKLNL MYIDSIDLEKITETEDPVKFHEAWQK LCKADGILVPGGFGIRGTLGKLQAISWARTKKIPFLGVCLGMQLAVIEFARNCLNLKDAD STEFRPNAPVPLVIDMPEHNPGNLGGTMRLGIRRTVFKTENSILRKLYGDVPFIEERHRH RFEVNPNIKQFEQNDLSFVGQDVGDRMEIIELANHPYFVG VQFHPEFSSRPMKPSPPY LGLLLAATGNLNAYLQQGCKLSSSDRYSDASDDSFSEPRIAELEIS


DKFZp686C17207, FLJ43358, MGC32997

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~65.70 kDa using CTPS2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CTPS2 expression in transfected 293T cell line by CTPS2 rabbit polyclonal antibody.

Lane 1: CTPS2 transfected lysate(65.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619