



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-CYP2D6-211

Anti CYP2D6

Product type	Primary antibodies
Immunogen	CYP2D6 (AAH75023.1, 1 a.a. ~ 497 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MGLEALVPLAVIVAIFLLLVDLMHRRQRWAARYPPGPLPLPGLGNLLHVDQNTPTCYCFDQ LRRRFGDVFSLQLAWTPVVVLNGLAAVREALVTHGEDTADRPPVPITQILGFGPRSQGVF LARYGPAWREQRRFSVSTLRNLGLGKKSLEQWVTEEAACLCAAFANHSGRPFRPNGLLDK AVSNVIASLTGRRFEYDDPRFLRLDLAQEGLKEESGFLREVLNAVVPVLLHIPALAGKV LRFQKAFLTQLDELLTEHRMTWDPAQPPRDLTEAFLAEMEKAKGNPESSFNDENLRIVVA DLFSAGMVTSTTLAWGLLLMILHPDVQRRVQQEIDDVIGQVRRPEMGDQAHMPYTTAVI HEVQRFQDIVPLGVTHMTSRDIEVQGFRIKGTTLITNLSSVLKDEAVWEKPFRRHPEHF LDAQGHFVKPEAFLPFSAGRRACLGEPLARMELFFFTSLLQHFSSVPTGQPRPSHHGV FAFLVSPSPYELCAVPR
Gene Alias	CPD6, CYP2D, CYP2D@, CYP2DL1, MGC120389, MGC120390, P450-DB1 P450C2D, P450DB1
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~55.80 kDa using CYP2D6 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CYP2D6 expression in transfected 293T cell line by CYP2D6 rabbit polyclonal antibody.

Lane 1: CYP2D6 transfected lysate(55.80 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619