

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CYP4X1-353

Anti CYP4X1

Product type	Primary antibodies
Immunogen	CYP4X1 (NP_828847.1, 1 a.a. ~ 509 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence

MEFSWLETRWARPFYLA FVFC LALG LLQA IKLYLRRQ RLLRDLRPFPAPPTHWFLGHQKF
IQDDNMEKLEEIIEKYPRAFPFWIGPFQAFFCIYDPDYAKTLLSRTDPKSQYLQKFSPPL
LGKGLAALDGPKEFQHRRLTPGFHFNILKAYIEVMAHSV KMMLDKWEKICSTQDTSVEV
YEHINSMSLDIIMKCAFSKETNCQTNSTHDPYAKAIFELSKIIFHRLYSLLYHSDIIFKL
SPQGYRFQKLSRVLNQYTDIIQERKKS LQAGVKQDNTPKRKYQDFLDIVLSAKDESGSS
FSDIDVHSEVSTFLLAGHDTLAASISWILYCLALNPEHQERCREEVRGILGDGSSITWDQ
LGEMSYTTMCIKETCRLIPAVPSISRDL SKPLTFPDGCTLPAGITVVL SIWGLHHNPAVW
KNPKVFDPLRFSQENS DQRHPYAYLPFSAGSRNCIGQEFAMIELKVTIALILLHFRVTPD
PTRPLTFPNHFILKPKNGMYLHLKKLSEC


MGC40051

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~58.90 kDa using CYP4X1 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of CYP4X1 expression in transfected 293T cell line by CYP4X1 rabbit polyclonal antibody.

Lane 1: CYP4X1 transfected lysate(58.90 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of CYP4X1 expression in human kidney (Predicted Mw: 58.3500 kDa).

Western Blotting (Tissue lysate)



CYP4X1 rabbit polyclonal antibody. Western Blot analysis of CYP4X1 expression in human kidney.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619