

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP-DHDH-549

Anti DHDH

Product type	Primary antibodies
Immunogen	DHDH (NP_055290.1, 1 a.a. ~ 334 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MALRWGIVSVGLISSDFTAVLQTLPRSEHQVVAVAARDLSRAKEFAQKHDIPKAYGSYEE LAKDPSVEVAYIGTQHPQHKAAMVLCCLAGKAVLCEKPTGVNAAEVREMVAEARSRALFL MEAIWTRFFPASEALRSVLAQGTGLDLRVARAEFGKNLIHVPRAVDRAQAGGALLDIGIY CVQFTSMVFGGQKPEKISVVGRRHETGVDDTVTVLLQYPGEVHGSFTCSITVQLSNTASV SGTKGMVQLLNPWCPTLVVKGEHKEFPLPPVPKDCNFDNGAGMSYEAHVWECLRKGM KESPVIPLSESELLADILEEVRKAIGVTFPQDKR

Gene Alias	HUM2DD
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~36.40 kDa using DHDH transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of DHDH expression in transfected 293T cell line by DHDH rabbit polyclonal antibody.

Lane 1: DHDH transfected lysate(36.40 KDa).
Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of DHDH transfected lysate

Immunoprecipitation



Immunoprecipitation of DHDH transfected lysate using DHDH rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with DHDH purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619